

# SCLA (Skiin Connected Life App) - Complete Reconstruction Documentation

---

## Executive Summary

---

This document provides a comprehensive analysis and reconstruction guide for the SCLA (Skiin Connected Life App) mobile application based on detailed screen analysis. The app is a health monitoring platform that integrates with SKIIN wearable devices to collect ECG data, track symptoms, and provide clinical insights for Holter studies.

## Key Features

- **Real-time ECG monitoring** with dual-channel display
- **Symptom logging** with intensity tracking and trigger identification
- **Blood pressure recording** with dual measurement capability
- **Diary functionality** for historical data review
- **Device management** for SKIIN pod pairing and monitoring
- **Clinical integration** for Holter study progress tracking

## Table of Contents

---

1. [App Architecture Overview](#)
2. [Screen-by-Screen Analysis](#)
3. [Navigation Flows](#)
4. [Component Hierarchy](#)
5. [Data Models](#)
6. [API Specifications](#)


7. [Technical Implementation Guide](#)
8. [User Experience Guidelines](#)

## App Architecture Overview

---

The SCLA app follows a tab-based navigation pattern with three main sections:

### Primary Navigation Structure



### Core Modules

1. **Authentication Module:** User sign-in with multiple methods
2. **Device Management Module:** SKIIN pod pairing and monitoring
3. **Health Data Module:** ECG, symptoms, and blood pressure tracking
4. **Analytics Module:** Data analysis and clinical insights
5. **Settings Module:** User preferences and device configuration

# Screen-by-Screen Analysis

---

## Authentication Screens

### Welcome Screen

- **Purpose:** App introduction and entry point
- **Key Elements:**
  - SKIIN branding and feature overview
  - "Get Started" button for new users
  - "Sign In" link for existing users
- **Navigation:** Entry point to registration or sign-in flows

### Sign In Screen

- **Purpose:** User authentication with multiple options
- **Key Elements:**
  - Email/password fields with show/hide toggle
  - "Forgot Password" link
  - Alternative authentication methods (6-digit code, QR code)
  - "Sign Up" link for new users
- **Validation:** Email format validation, password requirements
- **Security:** Secure credential handling, session management

## Device Pairing Screens

### Device Pairing Introduction

- **Purpose:** Guide users through pod pairing process
- **Key Elements:**
  - Visual instructions with device imagery
  - Step-by-step pairing guidance
  - "Pair Pod" and "Cancel" actions

- **Hardware Integration:** Bluetooth Low Energy (BLE) connectivity

## Pod Selection Screen

- **Purpose:** Select specific device from available options
- **Key Elements:**
  - List of discoverable devices with serial numbers
  - Radio button selection interface
  - Connection strength indicators
- **Technical Requirements:** BLE device discovery and filtering

## Connection Status Screen

- **Purpose:** Provide feedback during pairing process
- **Key Elements:**
  - Progress indication
  - Device information display
  - Connection status updates
- **User Experience:** Clear feedback and error handling

## Main Application Screens

### Home Dashboard

- **Purpose:** Central hub for status overview and quick access
- **Key Components:**
  - Greeting message with time-based personalization
  - Garment Signal Status card with real-time updates
  - Holter Study Progress tracking
  - 14-Day Progress visualization
  - Battery optimization prompts
- **Real-time Features:** Live signal status, sync indicators
- **Responsive Design:** Adaptive layout for different screen sizes

## ECG Viewer (Signal Status Screen)

- **Purpose:** Real-time ECG monitoring and signal quality assessment
- **Technical Specifications:**
  - Dual-channel ECG display (Channel 1 and Channel 3)
  - Configurable scale settings (mm/mV)
  - Real-time waveform rendering at 250Hz+ sample rate
  - Signal quality indicators with color coding
  - Heart rate calculation and display
- **User Interactions:**
  - Scale adjustment dropdown
  - Signal troubleshooting guidance
  - Navigation back to dashboard
- **Performance Requirements:** Low-latency data streaming, smooth animations

## Log Tab Interface

- **Purpose:** Data entry hub for health metrics
- **Modal Design:** "What would you like to add?" selection
- **Entry Types:**
  - Symptom logging with detailed attributes
  - Blood pressure recording with dual measurements
- **User Experience:** Quick access, minimal friction data entry

## Symptom Logging Flow

### Basic Symptom Entry

- **Date/Time Selection:** Calendar and time picker components
- **Symptom Selection:** Grid layout with predefined options
- **Search Functionality:** Real-time filtering of symptom list
- **Custom Symptoms:** Ability to add user-defined symptoms

## Detailed Symptom Entry

- **Intensity Rating:** 0-10 slider with visual feedback
- **Trigger Selection:** Multi-select with predefined and custom options
- **Duration Tracking:** Ongoing vs. intermittent with time inputs
- **Notes Field:** Free-text additional details
- **Validation:** Required field checking, data format validation

## Symptom Confirmation

- **Success Feedback:** Visual confirmation with illustration
- **ECG Analysis Integration:** Automatic correlation with ECG data
- **Analysis Timeline:** Expected completion time communication
- **Return Navigation:** Clear path back to main interface

## Blood Pressure Logging Flow

### Measurement Entry

- **Dual Readings:** First and second measurement capability
- **Numeric Input:** Custom keypad for precise entry
- **Measurement Guidelines:** Visual instructions for proper technique
- **Data Validation:** Range checking, format validation

### Confirmation and Storage

- **Success Confirmation:** Visual feedback for completed entry
- **Data Integration:** Automatic sync with health record
- **Historical Tracking:** Integration with diary timeline

## Diary Interface

### Calendar View

- **Monthly Navigation:** Month/year selection dropdown

- **Date Selection:** Interactive calendar with entry indicators
- **Entry Timeline:** Chronological list of daily entries
- **Refresh Capability:** Pull-to-refresh for data updates

## Entry Display

- **Card-Based Layout:** Distinct cards for different data types
- **Timestamp Precision:** Exact time display for each entry
- **Status Indicators:** Analysis and review status badges
- **Quick Actions:** Edit and delete functionality

## Entry Detail View

- **Comprehensive Display:** All entry attributes visible
- **Edit Functionality:** In-place editing capability
- **Remove Options:** Confirmation-based deletion
- **Analysis Status:** Pending/completed analysis indicators

## Settings and Management

### Device Management

- **Connection Status:** Real-time device connectivity
- **Sync Information:** Last sync time and current status
- **Battery Monitoring:** Device charge level display
- **Device Details:** Firmware, hardware version information
- **Unpair Functionality:** Device removal capability

### Advanced Settings

- **User Profile:** Personal information management
- **Clinical Program:** Study-specific configurations
- **Data Upload Preferences:** WiFi vs. cellular options
- **Custom Content Management:** Symptoms and triggers

## Manage Symptoms & Triggers

- **Custom Symptoms:** User-created symptom list
- **Custom Triggers:** User-defined trigger management
- **Delete Functionality:** Remove custom entries
- **Sync Integration:** Cloud synchronization of custom data

## Navigation Flows

---

### Primary User Journeys

#### First-Time User Flow

Welcome Screen → Sign Up → Device Pairing → Tutorial → Home Dashboard

#### Returning User Flow

Sign In → Home Dashboard → [Daily Usage Patterns]

#### Symptom Logging Journey

Home/Log Tab → Add Data Modal → Symptom Entry → Details Form → Confirmation → Diary Update

#### ECG Monitoring Journey

Home Dashboard → Signal Status → ECG Viewer → Troubleshooting (if needed)

#### Data Review Journey

Diary Tab → Calendar Navigation → Entry Selection → Detail View → Edit/Remove



## Modal and Popup Interactions

### Battery Optimization Popup

- **Trigger:** System detection of suboptimal settings
- **Options:** Accept optimization or dismiss
- **Integration:** Direct link to device settings

### Signal Quality Alerts

- **Trigger:** Poor ECG signal detection
- **Actions:** Troubleshooting guidance or dismissal
- **Context:** Real-time signal monitoring

### Add Data Modal

- **Trigger:** Log tab access or plus button
- **Options:** Symptom or blood pressure entry
- **Design:** Bottom sheet modal with clear options

## Component Hierarchy

---

### UI Component Structure

#### Navigation Components

- **BottomTabNavigator:** Main app navigation
- **TopNavigationBar:** Context-sensitive header
- **SideNavigationMenu:** Settings and profile access
- **BackButton:** Consistent navigation pattern

#### Data Display Components

- **ECGViewer:** Real-time waveform display
- **SignalStatusIndicator:** Connection and quality status

- **EntryCard:** Diary entry display component
- **ProgressIndicator:** Study and sync progress
- **CalendarComponent:** Date selection and navigation

## Input Components

- **SymptomSelector:** Grid-based symptom selection
- **IntensitySlider:** 0-10 rating scale
- **TriggerSelector:** Multi-select trigger interface
- **DateTimePicker:** Date and time selection
- **NumericKeypad:** Custom numeric input
- **SearchInput:** Real-time filtering capability

## Modal Components

- **AddDataModal:** Entry type selection
- **ConfirmationModal:** Success feedback
- **AlertModal:** Error and warning messages
- **SettingsModal:** Configuration interfaces

## State Management Architecture

### Global State

- **User Authentication:** Login status, user profile
- **Device Connection:** Pairing status, sync state
- **Real-time Data:** ECG stream, signal quality
- **Application Settings:** Preferences, configurations

### Local State

- **Form Data:** Entry forms, validation states
- **UI State:** Modal visibility, navigation state
- **Cache Data:** Offline storage, sync queue

- **Temporary Data:** Search results, filtered lists

## Data Models

---

### Core Entity Relationships

```
User (1) ↔ (M) UserDevices ↔ (M) Devices
User (1) ↔ (M) Symptoms
User (1) ↔ (M) BloodPressureReadings
User (1) ↔ (M) ECGData
Symptoms (1) ↔ (M) SymptomTriggers
User (1) ↔ (M) CustomSymptoms
User (1) ↔ (M) CustomTriggers
User (1) ↔ (1) UserPreferences
User (1) ↔ (M) HolterStudies
ECGData (1) ↔ (M) ECGAnalyses
```

### Key Data Structures

#### User Profile

```
interface User {
  id: string;
  email: string;
  firstName: string;
  lastName: string;
  dateOfBirth: Date;
  createdAt: Date;
  lastLogin: Date;
  isActive: boolean;
}
```

#### Device Information

```
interface Device {
  id: string;
  serialNumber: string;
  deviceType: 'skin_pod' | 'chestband';
  firmwareVersion: string;
  hardwareVersion: string;
  batteryLevel: number; // 0-100
  connectionStatus: 'connected' | 'disconnected' | 'syncing';
  lastSyncAt: Date;
}
```

## Symptom Entry

```
interface Symptom {
  id: string;
  userId: string;
  symptomName: string;
  intensity: number; // 0-10
  experiencedAt: Date;
  durationType: 'ongoing' | 'intermittent';
  durationHours?: number;
  durationMinutes?: number;
  triggers: string[];
  notes?: string;
  analysisStatus: 'pending' | 'analyzing' | 'completed';
  reviewStatus: 'pending' | 'reviewed' | 'approved';
}
```

## ECG Data

```
interface ECGData {
  id: string;
  userId: string;
  deviceId: string;
  recordedAt: Date;
  durationSeconds: number;
  sampleRate: number; // Hz
  channel1Data: number[]; // Voltage values
  channel3Data: number[];
  signalQuality: {
    channel1: 'good' | 'poor' | 'no_signal';
    channel3: 'good' | 'poor' | 'no_signal';
  };
  heartRate?: number;
  analysisStatus: 'pending' | 'analyzing' | 'completed' | 'failed';
}
```

## Blood Pressure Reading

```
interface BloodPressureReading {
  id: string;
  userId: string;
  measuredAt: Date;
  systolic1: number;
  diastolic1: number;
  systolic2?: number;
  diastolic2?: number;
  notes?: string;
}
```

# API Specifications

---

## Authentication Endpoints

### POST /api/auth/login

#### Request:

```
{
  "email": "user@example.com",
  "password": "securePassword123"
}
```

#### Response:

```
{
  "accessToken": "eyJhbGciOiJIUzI1NiIs... ",
  "refreshToken": "eyJhbGciOiJIUzI1NiIs... ",
  "user": {
    "id": "550e8400-e29b-41d4-a716-446655440000",
    "email": "user@example.com",
    "firstName": "John",
    "lastName": "Doe"
  },
  "expiresIn": 900
}
```

### POST /api/auth/login/qr

#### Request:

```
{
  "qrCode": "data:image/png;base64,iVBORw0KGgoAAAANSUhEUgAA..."
}
```

### POST /api/auth/login/six-digit

#### Request:

```
{
  "code": "123456",
  "deviceId": "device_identifier"
}
```

## Device Management Endpoints

GET /api/devices/discover

Response:

```
{
  "devices": [
    {
      "serialNumber": "31067601890",
      "deviceType": "skiin_pod",
      "signalStrength": -45,
      "isAvailable": true
    },
    {
      "serialNumber": "31046900579",
      "deviceType": "skiin_pod",
      "signalStrength": -62,
      "isAvailable": true
    }
  ]
}
```

POST /api/devices/pair

Request:

```
{
  "serialNumber": "31067601890",
  "deviceName": "My Skiin Pod"
}
```

Response:

```
{
  "device": {
    "id": "550e8400-e29b-41d4-a716-446655440000",
    "serialNumber": "31067601890",
    "deviceName": "My Skiin Pod",
    "pairingStatus": "success",
    "firmwareVersion": "18.0.9.1",
    "hardwareVersion": "20"
  }
}
```

GET /api/devices/status

Response:

```
{
  "devices": [
    {
      "id": "550e8400-e29b-41d4-a716-446655440000",
      "serialNumber": "31067601890",
      "connectionStatus": "connected",
      "batteryLevel": 75,
      "lastSyncAt": "2025-06-24T11:12:00Z",
      "syncStatus": "completed",
      "signalQuality": {
        "channel1": "good",
        "channel3": "poor"
      }
    }
  ]
}
```

## Health Data Endpoints

GET /api/ecg/realtime

WebSocket Connection:

```
wss://api.scla.com/ws/ecg/{userId}
Authorization: Bearer {accessToken}
```

Data Stream Format:

```
{
  "timestamp": "2025-06-24T11:14:30.123Z",
  "deviceId": "550e8400-e29b-41d4-a716-446655440000",
  "sampleRate": 250,
  "channel1": [0.1, 0.15, 0.12, 0.08, ...],
  "channel3": [0.05, 0.08, 0.06, 0.04, ...],
  "heartRate": 72,
  "signalQuality": {
    "channel1": "good",
    "channel3": "poor"
  }
}
```

POST /api/symptoms

Request:

```
{
  "symptomName": "Brain Fog",
  "intensity": 5,
  "experiencedAt": "2025-06-24T11:13:00Z",
  "durationType": "ongoing",
  "durationHours": 2,
  "durationMinutes": 30,
  "triggers": ["Caffeine", "Not Sure", "hot"],
  "notes": "Felt foggy after morning coffee"
}
```

## Response:

```
{
  "symptom": {
    "id": "550e8400-e29b-41d4-a716-446655440000",
    "symptomName": "Brain Fog",
    "intensity": 5,
    "experiencedAt": "2025-06-24T11:13:00Z",
    "analysisStatus": "pending",
    "ecgCorrelationId": "550e8400-e29b-41d4-a716-446655440001"
  }
}
```

## GET /api/symptoms

**Query Parameters:** - `startDate`: ISO 8601 date string - `endDate`: ISO 8601 date string  
 - `page`: Page number (default: 1) - `limit`: Items per page (default: 20)

## Response:

```
{
  "symptoms": [
    {
      "id": "550e8400-e29b-41d4-a716-446655440000",
      "symptomName": "Brain Fog",
      "intensity": 5,
      "experiencedAt": "2025-06-24T11:13:00Z",
      "triggers": ["Caffeine", "Not Sure", "hot"],
      "analysisStatus": "completed",
      "reviewStatus": "pending"
    }
  ],
  "pagination": {
    "page": 1,
    "limit": 20,
    "total": 45,
    "totalPages": 3
  }
}
```



## POST /api/blood-pressure

### Request:

```
{
  "measuredAt": "2025-06-24T11:12:00Z",
  "systolic1": 128,
  "diastolic1": 84,
  "systolic2": 129,
  "diastolic2": 85,
  "notes": "Morning reading after coffee"
}
```

## Diary and Analytics Endpoints

### GET /api/diary

**Query Parameters:** - `date`: YYYY-MM-DD format - `month`: YYYY-MM format - `year`: YYYY format

### Response:

```
{
  "date": "2025-06-24",
  "entries": [
    {
      "id": "550e8400-e29b-41d4-a716-446655440000",
      "type": "symptom",
      "timestamp": "2025-06-24T10:38:00Z",
      "data": {
        "symptomName": "Sweating",
        "intensity": 5,
        "triggers": ["Caffeine"],
        "analysisStatus": "pending",
        "reviewStatus": "pending"
      }
    },
    {
      "id": "550e8400-e29b-41d4-a716-446655440001",
      "type": "blood_pressure",
      "timestamp": "2025-06-24T11:12:00Z",
      "data": {
        "systolic": 128.5,
        "diastolic": 84,
        "reading": "128.5/84 mmHg"
      }
    }
  ]
}
```

## GET /api/holter-study/progress

### Response:

```
{
  "study": {
    "id": "550e8400-e29b-41d4-a716-446655440000",
    "name": "14-Day Holter Study",
    "startDate": "2025-06-10",
    "endDate": "2025-06-24",
    "progressPercentage": 65,
    "daysCompleted": 9,
    "totalDays": 14,
    "status": "active",
    "lastUpdated": "2025-06-24T10:21:00Z"
  }
}
```

## Settings and Preferences Endpoints

### GET /api/user/preferences

### Response:

```
{
  "preferences": {
    "dataUploadPreference": "wifi_cellular",
    "notificationEnabled": true,
    "ecgScaleSetting": 10.0,
    "timezone": "America/New_York",
    "language": "en"
  }
}
```

### PUT /api/user/preferences

### Request:

```
{
  "dataUploadPreference": "wifi_only",
  "notificationEnabled": false,
  "ecgScaleSetting": 5.0
}
```

### GET /api/custom-symptoms

### Response:

```
{
  "customSymptoms": [
    {
      "id": "550e8400-e29b-41d4-a716-446655440000",
      "name": "Custom Headache Type",
      "createdAt": "2025-06-20T14:30:00Z"
    }
  ]
}
```

**POST /api/custom-triggers**

**Request:**

```
{
  "triggerName": "hot weather"
}
```

**DELETE /api/custom-triggers/{triggerId}**

## Technical Implementation Guide

---

### Mobile App Architecture

#### Technology Stack Recommendations

- **Framework:** React Native or Flutter for cross-platform development
- **State Management:** Redux Toolkit or MobX for complex state
- **Navigation:** React Navigation or Flutter Navigator
- **Real-time Communication:** WebSocket with Socket.io
- **Local Storage:** SQLite with encryption
- **Charts/Graphs:** Victory Native or FL Chart
- **Bluetooth:** React Native BLE Manager or Flutter Blue

## Project Structure

```
src/
├── components/           # Reusable UI components
│   ├── common/          # Generic components
│   ├── ecg/             # ECG-specific components
│   ├── forms/           # Form components
│   └── navigation/      # Navigation components
├── screens/             # Screen components
│   ├── auth/            # Authentication screens
│   ├── home/            # Home dashboard
│   ├── log/             # Data entry screens
│   ├── diary/           # Diary and history
│   └── settings/        # Settings and preferences
├── services/            # API and business logic
│   ├── api/             # API client and endpoints
│   ├── bluetooth/       # Device communication
│   ├── storage/         # Local data management
│   └── analytics/       # Data analysis services
├── store/               # State management
│   ├── slices/          # Redux slices or MobX stores
│   └── middleware/      # Custom middleware
├── utils/              # Utility functions
│   ├── validation/      # Form validation
│   ├── formatting/      # Data formatting
│   └── constants/       # App constants
└── assets/              # Images, fonts, etc.
```

## Key Implementation Considerations

### Real-time ECG Display

```
// ECG Viewer Component Implementation
interface ECGViewerProps {
  deviceId: string;
  scaleSetting: number;
  onSignalQualityChange: (quality: SignalQuality) => void;
}

const ECGViewer: React.FC<ECGViewerProps> = ({
  deviceId,
  scaleSetting,
  onSignalQualityChange
}) => {
  const [ecgData, setEcgData] = useState<ECGDataPoint[]>([]);
  const [signalQuality, setSignalQuality] = useState<SignalQuality>();

  useEffect(() => {
    const websocket = new WebSocket(`wss://api.scla.com/ws/ecg/${deviceId}`);

    websocket.onmessage = (event) => {
      const data = JSON.parse(event.data);
      setEcgData(prev => [...prev.slice(-1000), ...data.samples]);
      setSignalQuality(data.signalQuality);
      onSignalQualityChange(data.signalQuality);
    };

    return () => websocket.close();
  }, [deviceId]);

  return (
    <View style={styles.container}>
      <ECGChart
        data={ecgData}
        scale={scaleSetting}
        signalQuality={signalQuality}
      />
      <SignalQualityIndicator quality={signalQuality} />
    </View>
  );
};
```

## Bluetooth Device Management

```
// Device Service Implementation
class DeviceService {
  private bleManager: BleManager;

  async discoverDevices(): Promise<Device[]> {
    const devices = await this.bleManager.startDeviceScan(
      ['SKIIN_SERVICE_UUID'],
      { allowDuplicates: false }
    );

    return devices.filter(device =>
      device.name?.includes('SKIIN') &&
      device.isConnectable
    );
  }

  async pairDevice(serialNumber: string): Promise<PairingResult> {
    try {
      const device = await this.bleManager.connectToDevice(serialNumber);
      await device.discoverAllServicesAndCharacteristics();

      // Setup ECG data streaming
      await this.setupECGStreaming(device);

      return { success: true, device };
    } catch (error) {
      return { success: false, error: error.message };
    }
  }

  private async setupECGStreaming(device: Device): Promise<void> {
    await device.monitorCharacteristicForService(
      'ECG_SERVICE_UUID',
      'ECG_DATA_CHARACTERISTIC_UUID',
      (error, characteristic) => {
        if (characteristic?.value) {
          const ecgData = this.parseECGData(characteristic.value);
          this.onECGDataReceived(ecgData);
        }
      }
    );
  }
}
```

## Form Validation and State Management

```
// Symptom Form Implementation
interface SymptomFormState {
  symptomName: string;
  intensity: number;
  experiencedAt: Date;
  triggers: string[];
  durationType: 'ongoing' | 'intermittent';
  durationHours?: number;
  durationMinutes?: number;
  notes: string;
}

const useSymptomForm = () => {
  const [formState, setFormState] = useState<SymptomFormState>({
    symptomName: '',
    intensity: 0,
    experiencedAt: new Date(),
    triggers: [],
    durationType: 'ongoing',
    notes: ''
  });

  const [errors, setErrors] = useState<Record<string, string>>({});

  const validateForm = (): boolean => {
    const newErrors: Record<string, string> = {};

    if (!formState.symptomName.trim()) {
      newErrors.symptomName = 'Symptom name is required';
    }

    if (formState.intensity < 0 || formState.intensity > 10) {
      newErrors.intensity = 'Intensity must be between 0 and 10';
    }

    if (formState.durationType === 'intermittent') {
      if (!formState.durationHours && !formState.durationMinutes) {
        newErrors.duration = 'Duration is required for intermittent symptoms';
      }
    }

    setErrors(newErrors);
    return Object.keys(newErrors).length === 0;
  };

  const submitForm = async (): Promise<boolean> => {
    if (!validateForm()) return false;

    try {
      await apiClient.post('/api/symptoms', formState);
      return true;
    } catch (error) {
      setErrors({ submit: 'Failed to save symptom' });
      return false;
    }
  };

  return {
    formState,
  }
}
```

```
    setFormState,  
    errors,  
    validateForm,  
    submitForm  
  };  
};
```

## Backend Implementation

### Technology Stack Recommendations

- **Runtime:** Node.js with TypeScript or Python with FastAPI
- **Database:** PostgreSQL with TimescaleDB for time-series data
- **Real-time:** Socket.io or WebSocket with Redis for scaling
- **Authentication:** JWT with refresh tokens
- **File Storage:** AWS S3 or Google Cloud Storage for ECG data
- **Message Queue:** Redis or RabbitMQ for analysis jobs
- **Monitoring:** Prometheus with Grafana
- **Deployment:** Docker with Kubernetes

### Database Optimization

#### Time-Series Data Handling

```
-- TimescaleDB hypertable for ECG data  
CREATE TABLE ecg_data (  
  id UUID DEFAULT gen_random_uuid(),  
  user_id UUID NOT NULL,  
  device_id UUID NOT NULL,  
  recorded_at TIMESTAMPTZ NOT NULL,  
  channel_1_data BYTEA,  
  channel_3_data BYTEA,  
  heart_rate INTEGER,  
  signal_quality JSONB,  
  PRIMARY KEY (id, recorded_at)  
);  
  
-- Convert to hypertable for time-series optimization  
SELECT create_hypertable('ecg_data', 'recorded_at');  
  
-- Create indexes for common queries  
CREATE INDEX idx_ecg_user_time ON ecg_data (user_id, recorded_at DESC);  
CREATE INDEX idx_ecg_device_time ON ecg_data (device_id, recorded_at DESC);
```



## Data Retention Policies

```
-- Automatic data retention for ECG data
SELECT add_retention_policy('ecg_data', INTERVAL '2 years');

-- Compression for older data
SELECT add_compression_policy('ecg_data', INTERVAL '30 days');
```

## API Implementation Examples

### Real-time ECG Streaming

```
// WebSocket handler for ECG streaming
class ECGStreamHandler {
  private connectedClients = new Map<string, WebSocket>();

  handleConnection(ws: WebSocket, userId: string) {
    this.connectedClients.set(userId, ws);

    ws.on('close', () => {
      this.connectedClients.delete(userId);
    });

    // Send initial signal status
    this.sendSignalStatus(userId);
  }

  async broadcastECGData(userId: string, ecgData: ECGDataPoint) {
    const client = this.connectedClients.get(userId);
    if (client && client.readyState === WebSocket.OPEN) {
      client.send(JSON.stringify({
        type: 'ecg_data',
        data: ecgData,
        timestamp: new Date().toISOString()
      }));
    }
  }

  async sendSignalStatus(userId: string) {
    const status = await this.getLatestSignalStatus(userId);
    const client = this.connectedClients.get(userId);

    if (client && client.readyState === WebSocket.OPEN) {
      client.send(JSON.stringify({
        type: 'signal_status',
        data: status
      }));
    }
  }
}
```

## Symptom Analysis Service

```
// Symptom analysis and ECG correlation
class SymptomAnalysisService {
  async analyzeSymptom(symptomId: string): Promise<AnalysisResult> {
    const symptom = await this.getSymptom(symptomId);
    const ecgData = await this.getECGDataAroundTime(
      symptom.userId,
      symptom.experiencedAt,
      { beforeMinutes: 30, afterMinutes: 30 }
    );

    if (!ecgData.length) {
      return { status: 'no_ecg_data', confidence: 0 };
    }

    // Perform analysis
    const analysis = await this.performECGAnalysis(ecgData, symptom);

    // Store results
    await this.storeAnalysisResults(symptomId, analysis);

    // Notify user if significant findings
    if (analysis.confidence > 0.8) {
      await this.notifyUser(symptom.userId, analysis);
    }

    return analysis;
  }

  private async performECGAnalysis(
    ecgData: ECGDataPoint[],
    symptom: Symptom
  ): Promise<AnalysisResult> {
    // Implement ECG analysis algorithms
    // This would typically involve:
    // 1. Heart rate variability analysis
    // 2. Arrhythmia detection
    // 3. QT interval measurement
    // 4. Correlation with symptom timing

    return {
      status: 'completed',
      confidence: 0.85,
      findings: ['Normal sinus rhythm', 'No significant arrhythmias'],
      recommendations: ['Continue monitoring']
    };
  }
}
```

# Security Implementation

## Authentication and Authorization

```
// JWT token management
class AuthService {
  generateTokens(user: User): TokenPair {
    const accessToken = jwt.sign(
      { userId: user.id, email: user.email },
      process.env.JWT_SECRET!,
      { expiresIn: '15m' }
    );

    const refreshToken = jwt.sign(
      { userId: user.id, tokenType: 'refresh' },
      process.env.REFRESH_SECRET!,
      { expiresIn: '30d' }
    );

    return { accessToken, refreshToken };
  }

  async validateToken(token: string): Promise<User | null> {
    try {
      const payload = jwt.verify(token, process.env.JWT_SECRET!) as JWPayload;
      return await this.getUserById(payload.userId);
    } catch (error) {
      return null;
    }
  }
}
```

## Data Encryption

```
// Health data encryption
class EncryptionService {
  private algorithm = 'aes-256-gcm';

  encryptHealthData(data: any): EncryptedData {
    const key = crypto.scryptSync(process.env.ENCRIPTION_KEY!, 'salt', 32);
    const iv = crypto.randomBytes(16);
    const cipher = crypto.createCipher(this.algorithm, key, iv);

    let encrypted = cipher.update(JSON.stringify(data), 'utf8', 'hex');
    encrypted += cipher.final('hex');

    const authTag = cipher.getAuthTag();

    return {
      encrypted,
      iv: iv.toString('hex'),
      authTag: authTag.toString('hex')
    };
  }

  decryptHealthData(encryptedData: EncryptedData): any {
    const key = crypto.scryptSync(process.env.ENCRIPTION_KEY!, 'salt', 32);
    const decipher = crypto.createDecipher(
      this.algorithm,
      key,
      Buffer.from(encryptedData.iv, 'hex')
    );

    decipher.setAuthTag(Buffer.from(encryptedData.authTag, 'hex'));

    let decrypted = decipher.update(encryptedData.encrypted, 'hex', 'utf8');
    decrypted += decipher.final('utf8');

    return JSON.parse(decrypted);
  }
}
```

## User Experience Guidelines

---

### Design Principles

#### Visual Design

- **Color Scheme:**
- Primary: Blue (#1976D2) for navigation and primary actions
- Secondary: Green (#4CAF50) for positive status indicators
- Warning: Orange (#FF9800) for attention-required states

- **Error:** Red (#F44336) for error states and alerts
- **Background:** Light gray (#F5F5F5) for main background
- **Cards:** White (FFFFFF) with subtle shadows

## Typography

- **Headers:** Bold, 18-24px for screen titles
- **Body Text:** Regular, 14-16px for content
- **Captions:** Light, 12-14px for timestamps and metadata
- **Font Family:** System fonts (San Francisco on iOS, Roboto on Android)

## Spacing and Layout

- **Margins:** 16px standard margin for screen edges
- **Padding:** 12px for card content, 8px for compact elements
- **Grid System:** 8px base unit for consistent spacing
- **Touch Targets:** Minimum 44px for interactive elements

## Accessibility Guidelines

### Screen Reader Support

- **Semantic Labels:** All interactive elements have descriptive labels
- **Content Description:** Complex UI elements have detailed descriptions
- **Navigation Hints:** Clear indication of navigation structure
- **State Announcements:** Changes in app state are announced

### Visual Accessibility

- **Color Contrast:** WCAG AA compliance (4.5:1 ratio minimum)
- **Text Scaling:** Support for dynamic type sizing
- **Focus Indicators:** Clear visual focus for keyboard navigation
- **Alternative Text:** All images have descriptive alt text

## Motor Accessibility

- **Touch Target Size:** Minimum 44x44pt touch targets
- **Gesture Alternatives:** Alternative input methods for complex gestures
- **Timeout Extensions:** Configurable timeouts for timed interactions
- **Error Prevention:** Clear validation and confirmation dialogs

## Performance Guidelines

### Loading and Response Times

- **App Launch:** < 3 seconds to interactive state
- **Screen Transitions:** < 300ms animation duration
- **API Responses:** < 2 seconds for data loading
- **Real-time Updates:** < 100ms latency for ECG streaming

### Memory and Battery Optimization

- **Memory Usage:** < 100MB baseline memory footprint
- **Battery Impact:** Optimized Bluetooth scanning and data processing
- **Background Processing:** Minimal background activity
- **Data Caching:** Intelligent caching to reduce network requests

### Offline Functionality

- **Core Features:** Basic app functionality available offline
- **Data Sync:** Automatic sync when connection restored
- **Conflict Resolution:** Clear handling of data conflicts
- **Storage Management:** Automatic cleanup of old cached data

## Error Handling and Recovery

### Error States

- **Network Errors:** Clear messaging with retry options

- **Device Connection:** Step-by-step troubleshooting guidance
- **Data Validation:** Inline validation with helpful error messages
- **System Errors:** Graceful degradation with fallback options

## Recovery Mechanisms

- **Automatic Retry:** Intelligent retry logic for transient failures
- **Manual Recovery:** Clear recovery actions for user-initiated fixes
- **Data Recovery:** Backup and restore capabilities for critical data
- **Support Integration:** Easy access to help and support resources

## Conclusion

---

This comprehensive documentation provides a complete blueprint for reconstructing the SCLA (Skiin Connected Life App) mobile application. The analysis covers all aspects from user interface design and navigation flows to backend architecture and data models.

## Key Implementation Priorities

1. **Real-time ECG Monitoring:** Core functionality requiring robust WebSocket implementation and optimized data visualization
2. **Device Integration:** Reliable Bluetooth connectivity with comprehensive error handling
3. **Health Data Management:** Secure, HIPAA-compliant storage and processing of sensitive health information
4. **User Experience:** Intuitive interface design with accessibility and performance optimization
5. **Clinical Integration:** Seamless workflow for healthcare providers and clinical studies

## Next Steps for Implementation

1. **Technical Architecture Setup:** Establish development environment and core infrastructure

2. **Authentication System:** Implement secure user authentication with multiple sign-in methods
3. **Device Pairing Flow:** Develop Bluetooth connectivity and device management
4. **Core UI Components:** Build reusable components for consistent user experience
5. **Real-time Data Streaming:** Implement ECG data visualization and processing
6. **Health Data Entry:** Create symptom and blood pressure logging functionality
7. **Diary and Analytics:** Develop historical data viewing and analysis features
8. **Testing and Validation:** Comprehensive testing including clinical validation
9. **Deployment and Monitoring:** Production deployment with monitoring and analytics

This documentation serves as a complete reference for development teams to recreate the SCLA app with full functionality and clinical-grade reliability.