

MVCP (Myant Virtual Clinic Portal) - Complete Application Reconstruction

Overview

The Myant Virtual Clinic Portal (MVCP) is a comprehensive web application designed for clinicians to manage patients, monitor Holter studies, review ECG data, and track symptoms. This document provides a complete reconstruction of the application based on screenshot analysis.

Application Structure

The MVCP application consists of several main sections accessible through a navigation menu:

1. **Symptom Log** - Main dashboard for reviewing patient symptoms
2. **Patients** - Patient management and registration
3. **Studies** - Holter study management and monitoring
4. **ECG Viewer** - Advanced ECG data visualization and analysis
5. **Notes** - Clinical notes and documentation
6. **Organization** - Team and member management

Screen-by-Screen Analysis

1. ECG Viewer - Symptom Annotation Screen

File: MVCP2025-06-24at10.39.26.png

Layout:

[← Back]

Annotations **for** symptom log

For symptom log on 2025-06-24, 10:38

[Save and

```
close |
```

[illegible]

Frequency filter	Mains noise filter	Amplitude	Date	Time
[0.05-45 Hz ▼]	[Off ▼]	[10mm/mV ▼]	[2025-06-24 ▼]	[10:34]
[Update]				[↺] [▶]
				ECG
Hours				

Full View: 10:34:30 - 10:38:30
Channel
☒ 1
☐ 2
☐ 3
◀ Previous. 4m
Next. 4m ▶

Expanded View: 10:38:20 - 10:38:30
All channels

[ECG waveform grid with detailed view]

Right Panel:

[Symptom Details]
[Contextual Metrics]

Sweating

Severity

5/10

Duration

On Going

Possible Triggers

Caffeine

Patient's note

Interactive Elements: - Back button: Returns to previous screen - Save and Close button: Saves annotations and closes modal - Observation checkboxes: Multiple selection for ECG observations - Notes text area: Free text input (1000 char limit) - Interpretation text area: Free text input (1000 char limit) - Filter controls: Dropdown menus for ECG display settings - Update button: Applies filter changes - Channel selection: Radio buttons for ECG channels - Navigation arrows: Previous/Next 4-minute segments - ECG Hours link: Likely opens time-based navigation - Symptom Details/Contextual Metrics tabs: Switch between views

2. ECG Viewer - Symptom Logs Tab

File: MVCP2025-06-24at10.39.14.png

Layout:

[Profile Icon]

Yang, Simon

[Verified]

PATIENT PROFILE

STUDIES

ECG

VIEWER

NOTES

[Symptom logs]

[ECG tags]

Status (All) ▼

Date (All) ▼

[Reset]

[?]

Page 1 of 1

SYMPTOM DATE AND TIME

REPORTED DATE AND TIME

STATUS ▼

SYMPTOMS

EXPERIENCED

ANNOTATIONS

2025-06-24 10:38

2025-06-24 10:38

[Pending review ▼]

Sweating

Add annotations

Holter Dates

Frequency filter

Mains noise filter

Amplitude

Date

Time

[0.05-45 Hz ▼]

[Off ▼]

[10mm/mV ▼]

[2025-06-24 ▼]

[10:35]

[Update]

[↺]

[▶]

ECG

Hours

Full View: 10:35:30 - 10:39:30

Channel ● 1 ○ 2 ○ 3

◀ Previous. 4m

Next. 4m ▶

Expanded View: 10:35:30 - 10:35:40

All channels

[ECG waveform grid with detailed view - Ch1 labeled]

Right Panel:

[Symptom Details] [Contextual Metrics]

Click on a symptom log to view more details about what the patient experienced.

Interactive Elements: - Symptom logs/ECG tags tabs: Switch between different views - Status filter dropdown: Filter symptoms by status (All, Pending review, etc.) - Date filter dropdown: Filter by date range - Reset button: Clear all filters - Help icon (?): Provides contextual help - Status dropdown in table: Change individual symptom status - "Add annotations" link: Opens annotation modal for specific symptom - "Holter Dates" link: Likely shows study date information - All ECG viewer controls same as previous screen - Symptom log rows: Clickable to view details in right panel

3. Notes Section - New Note Modal

File: MVCP2025-06-24at10.38.06.png

Layout:

Notes

Date (All) ▼ [Reset]

[🔒] New Note]

DATE AND TIME ▲ | CREATED BY | TITLE | NOTE

New Note

DATE AND TIME

2025-06-24 10:38

[Insert Title]

[≡] [≡] [≡] [≡] [•] [•] [B] [I] [U] Normal ▼ [A]

Insert note

i

For record keeping purposes, you will be unable to edit or delete this note once it is created.

[Cancel] [Create]

Interactive Elements: - New Note button: Opens note creation modal - Date filter dropdown: Filter notes by date - Reset button: Clear filters - Title input field: Text input for note title - Rich text editor toolbar: Text formatting options (alignment, lists, bold, italic, underline) - Style dropdown: Text style selection - Note content area: Rich text editor for note body - Cancel button: Close modal without saving - Create button: Save new note

4. Notes Section - Empty State

File: MVCP2025-06-24at10.38.01.png

Layout:

The screenshot shows a web interface for a patient's notes. At the top, a header bar contains a profile icon, the patient's name 'Yang, Simon [Verified]', and navigation links: 'PATIENT PROFILE | STUDIES | ECG VIEWER | NOTES |'. Below the header, the main content area is titled 'Notes'. It features a date filter dropdown set to 'Date (All)' with a downward arrow, a '[Reset]' button, and a '+ New' button with a lock icon. Below these controls is a table header with columns: 'DATE AND TIME ▲', 'CREATED BY', 'TITLE', and 'NOTE'. The table body is empty, and a message 'No Notes found for the given filters!' is displayed in the center.

Interactive Elements: - Same header controls as new note modal - Empty state message when no notes match current filters

5. ECG Viewer - ECG Tags Tab (Empty State)

File: MVCP2025-06-24at10.37.54.png

Layout:

[Profile Icon]

Yang, Simon [Verified]

PATIENT PROFILE | STUDIES | ECG VIEWER | NOTES |

[Symptom logs]

[ECG tags]

Status (All) ▼

Date (All) ▼

[Reset]

[?]

No Data found!

Frequency filter

Mains noise filter

Amplitude

Date

Time

[0.05-45 Hz ▼]

[Off ▼]

[10mm/mV ▼]

[2025-06-24 ▼]

[10:34]

[Update] |

[↺] [▶]

ECG

Hours |

Full View: 10:34:32 - 10:38:32

Channel ● 1 ○ 2 ○ 3

◀ Previous. 4m

Next. 4m ▶ |

Expanded View: 10:34:32 - 10:34:42

All channels

[ECG waveform grid with detailed view - Ch1 labeled]

Right Panel:

[Symptom Details]

[Contextual Metrics]

Click on a symptom log to view more details about what the patient experienced.

6. Studies Section - Holter Study Setup Modal

File: MVCP2025-06-24at10.37.46.png

Layout:

Studies

[Set up Holter study] |

Holter studies

Page 1 of 1 |

START DATE ▼ | STUDY DURATION | ... | STUDY PROGRESS | EDF DOWNLOAD

2025-06-24 10:21 | 14-days | ... | Avg. 0 hrs of usable ECG/day | -

Set up Holter study

Select the duration of the Holter study

☐ 24-hour

☐ 48-hour

☐ 72-hour

☐ 7-day

☒ 14-day

Select study start date

Year

Month

[2025 ▼]

[June ▼]

SUN

MON

TUE

WED

THU

FRI

SAT

1

2

3

4

5

6

7

8

9

10

11

12

13

14

15

16

17

18

19

20

21

22

23

[24]

25

26

27

28

29

30

Select study start time

[10:37]

Holter Study setup information

Setup Type		
[In-Clinic Setup ▼]		
Garment Type	Garment Size	
[Skiip-Band ▼]	[XS ▼]	

Interactive Elements: - Set up Holter study button: Opens study creation modal - Duration radio buttons: Select study length (24-hour, 48-hour, 72-hour, 7-day, 14-day) - Year/Month dropdowns: Navigate calendar - Calendar grid: Select start date (24th is selected) - Time input: Set start time - Setup Type dropdown: Choose setup method - Garment Type dropdown: Select device type - Garment Size dropdown: Select size

7. Studies Section - Study Progress Modal

File: MVCP2025-06-24at10.37.37.png

Layout:

Studies

[Set up Holter study] |

Holter studies

1 of 1 |

Page

START DATE ▼	STUDY DURATION	...	STUDY PROGRESS	EDF DOWNLOAD
2025-06-24 10:21	14-days	...	Avg. 0 hrs of usable ECG/day	-

Study progress **for** 14-days Holter study started on 2025-06-24, 10:21 |

Select a date on the calendar to view an hour by hour breakdown of the patient's ECG data collection progress.

Usable ECG / Total ECG

Year

Month

[2025 ▼]

[June ▼]

SUN

MON

TUE

WED

THU

FRI

SAT

1

2

3

4

5

6

7

8

9

10

11

12

13

14

15

16

17

18

19

20

21

22

23

[24]

25

26

27

28

29

30

00:00

0 mins / 0 mins

12:00

0 mins / 0 mins

01:00

0 mins / 0 mins

13:00

0 mins / 0 mins

02:00

0 mins / 0 mins

14:00

0 mins / 0 mins

03:00

0 mins / 0 mins

15:00

0 mins / 0 mins

04:00	0 mins / 0 mins	16:00	0 mins / 0 mins
05:00	0 mins / 0 mins	17:00	0 mins / 0 mins
06:00	0 mins / 0 mins	18:00	0 mins / 0 mins
07:00	0 mins / 0 mins	19:00	0 mins / 0 mins
08:00	0 mins / 0 mins	20:00	0 mins / 0 mins
09:00	0 mins / 0 mins	21:00	0 mins / 0 mins
10:00	0 mins / 0 mins	22:00	0 mins / 0 mins
11:00	0 mins / 0 mins	23:00	0 mins / 0 mins
Legend (Percentage of Usable Data Collected for Each Hour):			
[Gray] 0%	[Light Green] 1-24%	[Green] 25-50%	[Dark Green] >50%
0 hr 0 min of usable ECG / 0 hr 0 min of total ECG for the day			
[Close]			

8. Studies Section - Main List View

File: MVCP2025-06-24at10.37.31.png

Layout:

[Profile Icon] Yang, Simon [Verified]		PATIENT PROFILE STUDIES ECG	
VIEWER NOTES			
Studies		[Set up Holter study]	
Holter studies		Page	
1 of 1			
START DATE ▼ STUDY DURATION STATUS ESTIMATED END DATE ACTUAL END DATE			
STUDY PROGRESS EDF DOWNLOAD			
2025-06-24 10:21 14-days Ongoing 2025-07-08 -- Avg. 0 hrs of usable ECG/day -			

Interactive Elements: - Study progress modal: - Year/Month dropdowns: Navigate calendar - Calendar dates: Click to view hourly breakdown - Close button: Close modal - Hourly data display: Shows usable vs total ECG minutes - Color-coded legend: Visual indicator for data quality - Studies list: - Column headers: Sortable (START DATE has sort indicator) - Study progress link: Opens progress modal - Set up Holter study button: Opens study creation modal

9. Patient Profile Section

File: MVCP2025-06-24at10.37.24.png

Layout:

[Profile Icon] Yang, Simon [Verified] PATIENT PROFILE | STUDIES | ECG
VIEWER | NOTES |

Patient Profile [Edit] [SCLA
Login Code] |

Country

[Switzerland ▼]

First name

Last name

Date of birth

Sex

[Simon]

[Yang]

[1989-07-10]

[Male ▼]

Email

Phone

[simon.yang.ch@gmail.com]

[ +41 78 795 00 09]

Health insurance policy

Insurance number

[Avenir (Groupe Mutuel) ▼]

[7561234567897]

Medical record number

[]

Street address

City

Postal code

[Ackersteinstrasse 76] [Zurich]

[8049]

Assigned physician

[Peter Wood ▼]

[Archive Patient]
history) |

Associated Pod: 31067601890 (view

10. Patients Section - Create SCLA Account Modal (Form)

File: MVCP2025-06-24at10.37.13.png

Layout:

Approve registrations and manage patients under your organization's care.

Create SCLA Account

Create a Skin Connected Life App (SCLA) account **for** a patient
at your clinic.

* Country

[Canada ▼]

* First name * Last name * Date of birth Sex

[Enter patient's] [Enter patient's] [yyyy-mm-dd] [Unspecified ▼]

first name last name

* Email (SCLA account will be created under this email) *

[Enter patient's email address]

* Phone

[ +1]

* Health insurance policy * Health insurance policy number

[OHIP ▼] [Enter policy number]

* Version code * Expiry Date

[Enter version code] [yyyy-mm-dd]

Medical record number

[Enter patient's medical record number]			
Street address	City	Province	Postal code
[Enter street]	[Enter city]	[Select province ▼]	[Enter postal]
address			code]
Assigned physician			
[Select ▼]			
[Cancel] [Create Account]			

Interactive Elements: - Patient Profile: - Edit button: Enable form editing - SCLA Login Code button: Generate/view login code - All form fields: Editable when in edit mode - Archive Patient button: Archive patient record - View history link: Show pod association history - Create SCLA Account Modal: - All form fields with validation (asterisk indicates required) - Country dropdown: Select patient's country - Date inputs: Date picker for birth date and expiry - Phone field: Country code selector with phone input - Cancel button: Close modal without saving - Create Account button: Submit form and create account

11. Patients Section - Create SCLA Account Confirmation

File: MVCP2025-06-24at10.37.07.png

Layout:

[Patients Account] [Create SCLA Account]

Approve registrations and manage patients under your organization's care.

[Search of 36] [All Patients ▼] Page 1

PATIENT NAME ▲	DATE ADDED ▼	STATUS ▲	DATE OF BIRTH	HEALTHCARD NO.	VERSION
Yang, Simon	2025-06-24, 10:16	Verified	1989-07-10 (35)	7561234567897	--
CA-Andriod-One, Fozia	2025-06-16, 11:00	Verified	2000-01-20 (25)	1235567880	bb

[Additional patient rows...]

Create SCLA Account

Create a Skin Connected Life App (SCLA) account for a patient at your clinic. Ensure you have their consent first.

☒ The patient has provided consent to create an SCLA account on their behalf.

☒ I, Simon Yang, am taking responsibility for creating this patient's SCLA account.

[Cancel] [Continue]

12. Patients Section - Main List View

File: MVCP2025-06-24at10.37.00.png

Layout:

[Profile Icon] Simon Yang


Staging Clinic

373535


Turn on multi-factor authentication

≡ Symptom Log

 Patients

 Studies

 Organization

 User Manual

 Logout

Patients
Account] |

[Create SCLA

Approve registrations and manage patients under your organization's care.

[Search
of 36 |

] [All Patients ▼]

Page 1

PATIENT NAME ▲	DATE ADDED ▼	STATUS ▲	DATE OF BIRTH	HEALTHCARD NO. VERSION
Yang, Simon	2025-06-24, 10:16	Verified	1989-07-10 (35)	7561234567897 --
CA-Andriod-One, Fozia	2025-06-16, 11:00	Verified	2000-01-20 (25)	1235567880 bb
CA-Andriod, Fozia	2025-06-16, 10:27	Verified	2000-01-20 (25)	1235567880 bb
Prabhakaran, Praseena	2025-06-12, 16:06	Verified	1990-10-04 (34)	2433243434 ff
Rhino, Rick-HolterCloseTestOne	2025-06-11, 11:49	Verified	1985-06-11 (40)	9876543211 RR
Prabhakaran, Praseena	2025-06-10, 13:57	Verified	1989-10-16 (35)	1234567890 PP
MHDS 2.3.1 test, Praseena	2025-05-12, 15:10	Verified	1990-10-04 (34)	1234567890 PP
android, Fozia 2.3.1	2025-05-12, 12:29	Verified	2000-01-20 (25)	1236567880 bb
MHDS 2.3 Study Test, Praseena	2025-04-23, 14:25	Verified	1001-08-08 (1023)	1234567890 zz
MHDS 2.3 Regression Test, Praseena	2025-04-23, 13:01	Verified	1989-10-16 (35)	1234567890 VV

Interactive Elements: - Left sidebar navigation: - Symptom Log: Navigate to symptom log section - Patients: Current section (highlighted) - Studies: Navigate to studies section - Organization: Navigate to organization management - User Manual: Open help documentation - Logout: Sign out of application - Patient list: - Search field: Filter patients by name - All Patients dropdown: Filter by patient status - Column headers: Sortable (arrows indicate sort direction) - Patient rows: Clickable to view patient details - Pagination: Navigate through patient pages - Create SCLA Account confirmation: - Consent checkboxes: Required confirmations - Cancel/Continue buttons: Modal navigation

13. Symptom Log - Main View with Status Filter

File: MVCP2025-06-24at10.36.51.png

Layout:

Symptom Log

Symptom logs reported **in** the last 6 months.

Status (All) ▲ Date (All) ▼ ☐ Only show logs with Routine ECG
[Reset] |

Page 1

of 21 |
| ☒ Pending review |
| ☒ MD to review |
| ☒ Reviewed by MD |
| ☒ Reviewed by tech |
| [Cancel] [Apply] |

PATIENT NAME	SYMPTOM DATE	SYMPTOMS EXPERIENCED	PHYSICIAN	TAGS
MHDS 2.3.1 test, Praseena	2025-05-13 12:54	testing	① Not Assigned	--
[Pending review]				
MHDS 2.3 Study Test, Praseena	2025-04-23 15:00	study started	Prasi MD	test --
[Pending review]				
MHDS 2.3 Regression Test, Praseena	2025-04-23 14:00	new patient new study	① Not Assigned	--
[Pending review]				
MHDS 2.3 Testing, Praseena	2025-04-23 12:59	testing	① Not Assigned	-
[Pending review]				
MHDS 2.3 test, Praseena	2025-04-22 23:25	Routine ECG	① Not Assigned	--
[Pending review]				
Android 1, Fozia 2.3	2025-04-21 12:25	Routine ECG	① Not Assigned	--
[Pending review]				
Android 1, Fozia 2.3	2025-04-21 12:23	band to left	① Not Assigned	--
[Pending review]				
Android 1, Fozia 2.3	2025-04-21 12:21	band normal	① Not Assigned	--
[Pending review]				
Android 1, Fozia 2.3	2025-04-21 12:14	moving band to right	① Not Assigned	--
[Pending review]				

14. Symptom Log - Log History Modal

File: MVCP2025-06-24at10.36.41.png

Layout:

Symptom Log

Symptom logs reported **in** the last 6 months.

Status (All) ▼ Date (All) ▼ ☐ Only show logs with Routine ECG
[Reset] |

Page 1
of 21 |

PATIENT NAME	SYMPTOM DATE	SYMPTOMS EXPERIENCED	PHYSICIAN	TAGS	STATUS ▼
MHDS 2.3.1 test, Praseena	2025-05-13 12:54	testing	ⓘ Not Assigned	--	[Pending review]
					[Log History]
MHDS 2.3 Study Test, Praseena	2025-04-23 15:00	study started	Prasi MD	--	[Pending review]

Log History

[Most Recent ▼]

○ 2025-05-13, 12:54

Symptom Log created by Praseena MHDS 2.3.1 test

[Close]

15. Symptom Log - Main View (Sidebar Visible)

File: MVCP2025-06-24at10.36.25.png






Layout:

[Profile Icon] Simon Yang

Staging Clinic

373535

Turn on multi-factor authentication

- ≡ Symptom Log
-  Patients
-  Studies
-  Organization
-  User Manual
-  Logout

Symptom Log









Symptom logs reported in the last 6 months.

Status (All) ▼ Date (All) ▼ ☐ Only show logs with Routine ECG

[Reset] |

Page 1

of 21 |

PATIENT NAME	SYMPTOM DATE	SYMPTOMS EXPERIENCED	PHYSICIAN	TAGS
MHDS 2.3.1 test, Praseena	2025-05-13 12:54	testing	 Not Assigned	--
[Pending review]				
MHDS 2.3 Study Test, Praseena	2025-04-23 15:00	study started	Prasi MD	
test -- [Pending review]				
MHDS 2.3 Regression Test, Praseena	2025-04-23 14:00	new patient new study	 Not Assigned	--
[Pending review]				
MHDS 2.3 Testing, Praseena	2025-04-23 12:59	testing	 Not Assigned	-
[Pending review]				
MHDS 2.3 test, Praseena	2025-04-22 23:25	Routine ECG	 Not Assigned	--
[Pending review]				
Android 1, Fozia 2.3	2025-04-21 12:25	Routine ECG	 Not Assigned	--
[Pending review]				
Android 1, Fozia 2.3	2025-04-21 12:23	band to left	 Not Assigned	--
[Pending review]				
Android 1, Fozia 2.3	2025-04-21 12:21	band normal	 Not Assigned	--
[Pending review]				
Android 1, Fozia 2.3	2025-04-21 12:14	moving band to right	 Not Assigned	--
[Pending review]				

16. Organization Management

File: MVCP2025-06-24at10.39.49.png

Layout:

Organization Member] |

[Add

Add and manage members **in** your organization.

[Search of 18 |

] [All Members ▼]

Page 1

MEMBER NAME	EMAIL	ROLE	DATE ADDED
(invited)	fozia.noor+da@myant.ca	Desk Admin	--
(invited)	@gmail.com	Desk Admin	--
(invited)	gurpreet_test351@mailinator.com	Desk Admin	--
(invited)	gurpreet_test153@dispostable.com	Cardiac Tech	--
(invited)	gurpreet_test210@mailinator.com	MD	--
(invited)	pranali.patel+stagmd@myant.ca	MD	--
(invited)	augtestingmhp@dispostable.com	Cardiac Tech	--
(invited)	tetiana.kokorovets+ref2@myant.ca	Desk Admin	--
(invited)	ppatel@dispostable.com	Cardiac Tech	--

Interactive Elements: - Symptom Log: - Status filter dropdown with multi-select checkboxes - Date filter dropdown - Routine ECG checkbox filter - Reset button: Clear all filters - Status dropdown in each row: Change symptom status - Log History link: Opens history modal - Patient name links: Navigate to patient details - Organization: - Add Member button: Invite new organization members - Search field: Filter members by name/email - All Members dropdown: Filter by member type/status - Member rows: Show invitation status and roles

17. ECG Viewer - Contextual Metrics View

File: MVCP2025-06-24at10.39.34.png


Layout:


Frequency filter	Mains noise filter	Amplitude	Date	Time
[0.05-45 Hz ▼]	[Off ▼]	[10mm/mV ▼]	[2025-06-24 ▼]	[10:34]
[Update]				[↺] [▶]
Hours				ECG
Full View: 10:34:30 - 10:38:30 Channel <input checked="" type="radio"/> 1 <input type="radio"/> 2 <input type="radio"/> 3 ◀ Previous. 4m				
Next. 4m ▶				


Expanded View: 10:38:20 - 10:38:30 All channels				
[ECG waveform grid - Ch1]				
[ECG waveform grid - Ch2 with glucose trend line]				
[ECG waveform grid - Ch3]				


Right Panel - Contextual Metrics:


[Symptom Details]	[Contextual Metrics]
-------------------	----------------------


 BLOOD PRESSURE	
-- mmHg (closest before)	
-- mmHg (closest after)	

 OXYGEN SATURATION	
-- mmHg (closest before)	
-- mmHg (closest after)	

 BODY WEIGHT	
-- kg (closest before)	
-- kg (closest after)	

 GLUCOSE	
-- mmol/L (closest before)	
-- mmol/L (closest after)	

 BODY TEMPERATURE	
-- °C	

 ACTIVITY	
-- steps	
-- (last posture)	

18. Studies Section - Overview with Multiple Studies

File: MVCP2025-06-24at10.40.08.png

Layout:

[Profile Icon] Simon Yang


Staging Clinic


373535


Turn on multi-factor authentication


≡ Symptom Log

 Patients

 Studies

 Organization

 User Manual

 Logout

Studies

Add and manage studies **in** your organization.

Holter studies [Status (All) ▼]
of 33 |

Page 1

PATIENT NAME	START DATE ▼	STUDY DURATION	STATUS	ESTIMATED END DATE	ACTUAL END DATE	STUDY PROGRESS	EDF DOWNLOAD
Fozia 2.3, Android 1	2028-04-18 09:50	14-days	Closed by Staff	2028-05-02	2028-04-18	Avg. 0 hrs of usable ECG/day	-
Praseena, MHDS 2.3 test	2026-04-01 16:50	48-hours	Cancelled	2026-04-03	2025-04-11	0 hrs of usable ECG data, 0 hrs of total ECG data	-
Fozia 2.3, Android 1	2026-04-01 09:57	14-days	Closed by Staff	2026-04-15	2026-04-01	Avg. 0 hrs of usable ECG/day	-
Simon, Yang	2025-06-24 10:21	14-days	Ongoing	2025-07-08	--	Avg. 0 hrs of usable ECG/day	-
Fozia, CA-Andriod-One	2025-06-16 11:05	24-hours	Closed by Staff	2025-06-17	2025-06-16	0 hrs of usable ECG data, 0 hrs of total ECG data	-

Interactive Elements: - ECG Viewer Contextual Metrics: - Symptom Details/Contextual Metrics tabs: Switch between views - Contextual metrics display: Shows physiological data closest to ECG timepoint - Multiple ECG channels with overlay data (glucose trend on Ch2) - All standard ECG viewer controls - Studies Overview: - Status filter dropdown: Filter studies by status - Column headers: Sortable (START DATE shows sort indicator) - Study progress links: Open detailed progress views - Patient name links: Navigate to

patient details - Multiple study statuses: Ongoing, Closed by Staff, Cancelled - EDF
Download column: Access to raw data files

Summary of Interactive Elements and Navigation

The MVCP application provides a comprehensive interface for managing cardiac monitoring studies with the following key interactive patterns:

1. **Tab-based navigation** within sections (Symptom logs/ECG tags, Symptom Details/Contextual Metrics)
2. **Modal dialogs** for complex forms (Create SCLA Account, Set up Holter study, New Note, Study Progress)
3. **Dropdown filters** with multi-select capabilities for data filtering
4. **Sortable table columns** with visual indicators
5. **Rich text editing** capabilities in notes section
6. **Advanced ECG viewer** with multiple channels, filtering, and time navigation
7. **Contextual data display** showing physiological metrics alongside ECG data
8. **Status management** with dropdown controls for workflow progression

Component Interactions and Navigation Flows

Main Navigation Structure

The MVCP application uses a sidebar navigation pattern with the following main sections:

```
Main Navigation:  
├─ Symptom Log (≡ icon)  
├─ Patients (👥 icon)  
├─ Studies (📋 icon)  
├─ Organization (🏢 icon)  
├─ User Manual (❓ icon)  
└─ Logout (🚪 icon)
```

Top-Level Navigation Bar

When viewing a specific patient, a secondary navigation appears:

Patient Context Navigation:

- PATIENT PROFILE
- STUDIES
- ECG VIEWER
- NOTES

Detailed Interaction Mappings

1. Symptom Log Section

Main View Interactions: - **Status Filter Dropdown** → Opens multi-select filter with options: - Pending review - MD to review - Reviewed by MD - Reviewed by tech - **Date Filter Dropdown** → Date range selection - **Routine ECG Checkbox** → Filters to show only routine ECG logs - **Reset Button** → Clears all applied filters - **Patient Name Links** → Navigate to Patient Profile for that patient - **Status Dropdown (per row)** → Change individual symptom status - **Log History Link** → Opens Log History modal - **Pagination Controls** → Navigate through symptom log pages

Log History Modal: - **Sort Dropdown** → Change chronological order (Most Recent, Oldest First, etc.) - **Close Button** → Close modal and return to main view

2. Patients Section

Main List View Interactions: - **Create SCLA Account Button** → Opens account creation workflow - **Search Field** → Real-time patient name filtering - **All Patients Dropdown** → Filter by patient status/type - **Column Headers** → Sort patients by different criteria - **Patient Name Links** → Navigate to Patient Profile - **Pagination Controls** → Navigate through patient pages

Create SCLA Account Workflow: 1. **Create SCLA Account Button** → Opens consent confirmation modal 2. **Consent Checkboxes** → Required confirmations 3. **Continue Button** → Opens patient information form 4. **Form Fields** → Patient data entry with validation 5. **Create Account Button** → Submit and create account 6. **Cancel Button** → Exit workflow at any step

3. Patient Profile Section

Profile View Interactions: - **Edit Button** → Enable form editing mode - **SCLA Login Code Button** → Generate/display patient login code - **Form Fields** → Editable when in edit mode - **Dropdown Selectors** → Country, insurance, physician selection - **Archive Patient Button** → Archive patient record - **View History Link** → Show pod association history

4. Studies Section

Main List View Interactions: - **Set up Holter study Button** → Opens study creation modal - **Status Filter Dropdown** → Filter studies by status - **Column Headers** → Sort studies by different criteria - **Study Progress Links** → Opens study progress modal - **Patient Name Links** → Navigate to Patient Profile - **EDF Download Links** → Download study data files

Set up Holter Study Modal: - **Duration Radio Buttons** → Select study length (24h, 48h, 72h, 7-day, 14-day) - **Year/Month Dropdowns** → Navigate calendar - **Calendar Grid** → Select start date - **Time Input** → Set start time - **Setup Type Dropdown** → Choose setup method - **Garment Type/Size Dropdowns** → Select device configuration

Study Progress Modal: - **Year/Month Dropdowns** → Navigate calendar - **Calendar Dates** → Click to view hourly breakdown - **Hourly Data Display** → Shows usable vs total ECG minutes - **Close Button** → Return to studies list

5. ECG Viewer Section

Main Viewer Interactions: - **Symptom logs/ECG tags Tabs** → Switch between data views - **Filter Controls:** - Frequency filter dropdown → Adjust ECG frequency filtering - Mains noise filter dropdown → Noise reduction settings - Amplitude dropdown → Adjust ECG amplitude scaling - Date picker → Select viewing date - Time input → Set specific time - Update button → Apply filter changes - **Channel Selection** → Radio buttons for ECG channels (1, 2, 3) - **Navigation Controls:** - Previous/Next 4m buttons → Time-based navigation - ECG Hours link → Time-based navigation interface - **Symptom Details/Contextual Metrics Tabs** → Switch right panel view

Symptom Logs Tab: - **Status/Date Filters** → Filter symptom data - **Reset Button** → Clear filters - **Help Icon (?)** → Contextual help - **Add annotations Link** → Opens

annotation modal - **Holter Dates Link** → Study date information - **Symptom Rows** → Click to view details in right panel

Annotation Modal: - **Back Button** → Return to previous view - **Observation Checkboxes** → Multiple selection for ECG findings - **Notes Text Area** → Free text input (1000 char limit) - **Interpretation Text Area** → Clinical interpretation (1000 char limit) - **Save and Close Button** → Save annotations and close modal

6. Notes Section

Main View Interactions: - **New Note Button** → Opens note creation modal - **Date Filter Dropdown** → Filter notes by date - **Reset Button** → Clear filters - **Column Headers** → Sort notes

New Note Modal: - **Title Input** → Note title - **Rich Text Toolbar** → Text formatting options - **Style Dropdown** → Text style selection - **Content Area** → Rich text editor - **Cancel Button** → Close without saving - **Create Button** → Save new note

7. Organization Section

Main View Interactions: - **Add Member Button** → Invite new organization members - **Search Field** → Filter members by name/email - **All Members Dropdown** → Filter by member type/status - **Member Rows** → View member details and status

Modal and Popup Behaviors

Modal Types: 1. **Form Modals** (Create SCLA Account, Set up Holter study, New Note) - Overlay main content - Require user action to close - Include Cancel/Submit buttons - Form validation before submission

1. **Information Modals** (Study Progress, Log History)
2. Display read-only information
3. Single Close button
4. Click outside to close (implied)
5. **Annotation Modal** (ECG Viewer)
6. Complex form with multiple sections
7. Save and Close functionality

8. Back button navigation

Popup Behaviors: 1. **Dropdown Filters** - Click to open/close - Multi-select capability (Status filters) - Apply/Cancel buttons for complex filters

1. **Status Dropdowns**
2. Inline editing capability
3. Immediate save on selection
4. Visual feedback on status change

ECG Viewer Advanced Functionality

Time Navigation: - 4-minute segments with Previous/Next navigation - Full view shows 4-minute window - Expanded view shows 10-second detailed window - ECG Hours link provides time-based navigation interface

Channel Management: - 3 ECG channels available - Radio button selection - Individual channel viewing - "All channels" view in expanded mode

Filter Controls: - Real-time ECG signal processing - Frequency filtering (0.05-45 Hz standard) - Mains noise filtering options - Amplitude scaling (10mm/mV standard) - Update button applies all filter changes

Contextual Data Integration: - Right panel shows physiological metrics - Blood pressure, oxygen saturation, body weight, glucose, temperature, activity - "Closest before/after" temporal correlation - Glucose trend overlay on ECG Channel 2

Visual App Structure and Relationship Mapping

Application Hierarchy

```
MVCP (Myant Virtual Clinic Portal)
├── Authentication Layer
│   ├── Login
│   ├── Multi-factor Authentication
│   └── Session Management
├── Main Application Shell
│   ├── User Profile Header
│   │   ├── Profile Icon
│   │   ├── User Name & Clinic
│   │   ├── User ID
│   │   └── MFA Toggle
│   ├── Primary Navigation Sidebar
│   │   ├── Symptom Log (≡)
│   │   ├── Patients (👥)
│   │   ├── Studies (📋)
│   │   ├── Organization (🏢)
│   │   ├── User Manual (❓)
│   │   └── Logout (🚪)
│   ├── Main Content Area
│   │   ├── Section-Specific Content
│   │   └── Modal Overlay Layer
├── Patient Context Shell (when patient selected)
│   ├── Patient Header
│   │   ├── Patient Name
│   │   ├── Verification Status
│   │   └── Patient Context Navigation
│   │       ├── PATIENT PROFILE
│   │       ├── STUDIES
│   │       ├── ECG VIEWER
│   │       └── NOTES
│   └── Patient-Specific Content Area
```

Section-Level Structure

1. Symptom Log Section

```
Symptom Log
├── Filter Controls
│   ├── Status Multi-Select Filter
│   ├── Date Range Filter
│   ├── Routine ECG Checkbox
│   └── Reset Button
├── Data Table
│   ├── Patient Name (→ Patient Profile)
│   ├── Symptom Date
│   ├── Symptoms Experienced
│   ├── Physician Assignment
│   ├── Tags
│   └── Status Dropdown
├── Pagination Controls
└── Modals
    ├── Log History Modal
    │   ├── Sort Controls
    │   ├── History Timeline
    │   └── Close Button
```

2. Patients Section

```
Patients
├── Action Bar
│   ├── Create SCLA Account Button
│   ├── Search Field
│   └── Patient Type Filter
├── Data Table
│   ├── Patient Name (→ Patient Profile)
│   ├── Date Added
│   ├── Status
│   ├── Date of Birth
│   ├── Healthcare Number
│   └── Version
├── Pagination Controls
└── Modals
    ├── Create SCLA Account Workflow
    │   ├── Consent Confirmation
    │   │   ├── Consent Checkboxes
    │   │   ├── Cancel Button
    │   │   └── Continue Button
    │   └── Patient Information Form
    │       ├── Personal Information Fields
    │       ├── Contact Information Fields
    │       ├── Insurance Information Fields
    │       ├── Address Information Fields
    │       ├── Medical Information Fields
    │       ├── Cancel Button
    │       └── Create Account Button
    └── Account Creation Confirmation
```

3. Patient Profile Section

```
Patient Profile
├─ Action Bar
│   └─ Edit Button
│   └─ SCLA Login Code Button
├─ Patient Information Form
│   └─ Country Selection
│   └─ Personal Information
│       └─ First Name
│       └─ Last Name
│       └─ Date of Birth
│       └─ Sex
│   └─ Contact Information
│       └─ Email
│       └─ Phone
│   └─ Insurance Information
│       └─ Health Insurance Policy
│       └─ Insurance Number
│   └─ Medical Information
│       └─ Medical Record Number
│   └─ Address Information
│       └─ Street Address
│       └─ City
│       └─ Postal Code
│   └─ Care Team
│       └─ Assigned Physician
├─ Patient Actions
│   └─ Archive Patient Button
└─ Associated Data
    └─ Pod Association (with history link)
```

4. Studies Section

```
Studies
├── Action Bar
│   ├── Set up Holter study Button
│   └── Status Filter
├── Data Table
│   ├── Patient Name (→ Patient Profile)
│   ├── Start Date
│   ├── Study Duration
│   ├── Status
│   ├── Estimated End Date
│   ├── Actual End Date
│   ├── Study Progress (→ Progress Modal)
│   └── EDF Download
├── Pagination Controls
└── Modals
    ├── Set up Holter Study
    │   ├── Duration Selection
    │   ├── Start Date Calendar
    │   ├── Start Time Input
    │   └── Study Configuration
    │       ├── Setup Type
    │       ├── Garment Type
    │       └── Garment Size
    └── Study Progress Modal
        ├── Calendar Navigation
        ├── Date Selection
        ├── Hourly Data Display
        ├── Quality Legend
        └── Close Button
```

5. ECG Viewer Section

```
ECG Viewer
├── Tab Navigation
│   ├── Symptom logs Tab
│   └── ECG tags Tab
├── Filter Controls (per tab)
│   ├── Status Filter
│   ├── Date Filter
│   ├── Reset Button
│   └── Help Icon
├── Data Table (tab-specific)
│   ├── Symptom Date/Time
│   ├── Reported Date/Time
│   ├── Status Dropdown
│   ├── Symptoms Experienced
│   └── Actions
│       ├── Add annotations (→ Annotation Modal)
│       └── Holter Dates
├── ECG Viewer Controls
│   ├── Filter Controls
│   │   ├── Frequency Filter
│   │   ├── Mains Noise Filter
│   │   ├── Amplitude Control
│   │   ├── Date Picker
│   │   ├── Time Input
│   │   └── Update Button
│   ├── Navigation Controls
│   │   ├── Channel Selection (1, 2, 3)
│   │   ├── Previous/Next 4m
│   │   └── ECG Hours Link
│   └── View Controls
│       ├── Full View (4-minute window)
│       └── Expanded View (10-second detail)
├── ECG Display Area
│   ├── Multi-channel ECG Waveforms
│   ├── Time Axis
│   └── Amplitude Scaling
├── Right Panel
│   ├── Tab Navigation
│   │   ├── Symptom Details Tab
│   │   └── Contextual Metrics Tab
│   └── Content Area
│       ├── Symptom Information
│       │   ├── Symptom Type
│       │   ├── Severity
│       │   ├── Duration
│       │   ├── Triggers
│       │   └── Patient Notes
│       └── Contextual Metrics
│           ├── Blood Pressure
│           ├── Oxygen Saturation
│           ├── Body Weight
│           ├── Glucose
│           ├── Body Temperature
│           └── Activity Data
└── Modals
    ├── Annotation Modal
    │   ├── Navigation
    │   │   └── Back Button
    │   └── Observation Section
```

- Medical Observation Checkboxes
- Notes Section
 - Free Text Area (1000 chars)
- Interpretation Section
 - Clinical Interpretation Area (1000 chars)
- Actions
 - Save and Close Button

6. Notes Section

- Notes
 - Action Bar
 - New Note Button
 - Date Filter
 - Reset Button
 - Data Table
 - Date and Time
 - Created By
 - Title
 - Note Content
 - Pagination Controls
 - Modals
 - New Note Modal
 - Date and Time (auto-populated)
 - Title Input
 - Rich Text Editor
 - Formatting Toolbar
 - Alignment Controls
 - List Controls
 - Text Formatting (B, I, U)
 - Style Dropdown
 - Text Color
 - Content Area
 - Information Notice
 - Actions
 - Cancel Button
 - Create Button

7. Organization Section

- Organization
 - Action Bar
 - Add Member Button
 - Search Field
 - Member Type Filter
 - Data Table
 - Member Name
 - Email
 - Role
 - Date Added
 - Pagination Controls

Cross-Section Relationships

Patient-Centric Navigation Flow

```
Patient Selection (from any section)
↓
Patient Context Activated
↓
Patient Context Navigation Available:
├─ PATIENT PROFILE → Patient Profile Section
├─ STUDIES → Studies filtered for this patient
├─ ECG VIEWER → ECG data for this patient
├─ NOTES → Notes for this patient
```

Data Relationship Flow

```
Patient Creation (Patients Section)
↓
Patient Profile (Patient Profile Section)
↓
Study Creation (Studies Section)
↓
ECG Data Collection (automatic)
↓
Symptom Logging (patient-initiated)
↓
Symptom Review (Symptom Log Section)
↓
ECG Analysis (ECG Viewer Section)
↓
Clinical Notes (Notes Section)
```

Modal Interaction Patterns

```
Primary Action Buttons → Form Modals
├─ Create SCLA Account → Multi-step form workflow
├─ Set up Holter study → Configuration form
├─ New Note → Rich text editor
├─ Add annotations → Clinical observation form

Information Links → Display Modals
├─ Study Progress → Calendar with hourly data
├─ Log History → Timeline view
├─ View History → Pod association history

Status Changes → Inline Dropdowns
├─ Symptom Status → Workflow progression
├─ Study Status → Study lifecycle management
```

Component Reusability Map

Shared Components

Data Tables

- └─ Sortable Headers
- └─ Pagination Controls
- └─ Search/Filter Controls
- └─ Action Buttons

Form Components

- └─ Dropdown Selectors
- └─ Date Pickers
- └─ Text Inputs
- └─ Rich Text Editors
- └─ Validation Messages

Navigation Components

- └─ Tab Navigation
- └─ Breadcrumb Navigation
- └─ Sidebar Navigation
- └─ Modal Navigation

Filter Components

- └─ Multi-select Dropdowns
- └─ Date Range Pickers
- └─ Checkbox Filters
- └─ Reset Controls

Specialized Components

ECG Viewer

- └─ Waveform Renderer
- └─ Channel Selector
- └─ Time Navigation
- └─ Filter Controls
- └─ Overlay Data Display

Study Progress

- └─ Calendar Grid
- └─ Quality Indicators
- └─ Hourly Data Display
- └─ Legend Components

Patient Context

- └─ Patient Header
- └─ Verification Badge
- └─ Context Navigation
- └─ Associated Data Display

User Flows and Workflows

Primary User Roles

The MVCP system supports multiple user roles with different access levels: - **Desk Admin:** Administrative functions, patient management - **MD (Medical Doctor):** Clinical review, diagnosis, patient care - **Cardiac Tech:** Technical ECG analysis, device management - **Organization Admin:** User management, system configuration

Core User Workflows

1. Create Patient Account Workflow

Trigger: Clinician needs to register a new patient for cardiac monitoring

Flow Steps:

1. Navigate **to** Patients Section
 - └ Click "**Patients**" in sidebar navigation
2. Initiate Account Creation
 - └ Click "**Create SCLA Account**" button
3. Consent Verification
 - └ Review consent requirements
 - └ Check "**Patient has provided consent**" checkbox
 - └ Check "**I am taking responsibility**" checkbox
 - └ Click "**Continue**" button
4. Patient Information Entry
 - └ Select Country (required)
 - └ Enter Personal Information
 - └ First Name (required)
 - └ Last Name (required)
 - └ Date of Birth (required)
 - └ Sex (optional)
 - └ Enter **Contact** Information
 - └ Email (required - becomes SCLA **login**)
 - └ Phone (required)
 - └ Enter Insurance Information
 - └ Health Insurance Policy (required)
 - └ Policy Number (required)
 - └ Version Code (required)
 - └ **Expiry** Date (required)
 - └ Enter Medical Information
 - └ Medical Record Number (optional)
 - └ Enter Address Information
 - └ Street Address (optional)
 - └ City (optional)
 - └ Province (optional)
 - └ **Postal** Code (optional)
 - └ Assign Physician (optional)
5. Account Creation
 - └ Click "**Create Account**" button
 - └ **System** **validates** all required fields
 - └ **System** creates SCLA account
 - └ Patient appears in patients **list** with "**Verified**" status
6. **Post**-Creation Actions
 - └ Generate SCLA **Login** Code (**if** needed)
 - └ Navigate **to** Patient Profile **for** additional configuration

Error Handling: - Form validation prevents submission with missing required fields -
Email validation ensures proper format - Phone number validation with country code -
Insurance policy validation

2. Create Holter Study Workflow

Trigger: Patient requires cardiac monitoring study

Prerequisites: Patient account must exist and be verified

Flow Steps:

1. Navigate **to** Studies Section
 - └ Click "**Studies**" in sidebar navigation
2. Initiate Study Creation
 - └ Click "**Set up Holter study**" button
3. Study Configuration
 - └ Select Study Duration (required)
 - └ 24-hour
 - └ 48-hour
 - └ 72-hour
 - └ 7-day
 - └ 14-day (**default** shown)
 - └ Select Start Date (required)
 - └ Navigate **to** desired month/year
 - └ Click **on** calendar date
 - └ Set Start Time (required)
 - └ Enter time in HH:MM **format**
 - └ Configure Study Setup
 - └ Setup Type (required)
 - └ In-Clinic Setup (**default**)
 - └ Garment Type (required)
 - └ Skiip-Band (**default**)
 - └ Garment Size (required)
 - └ XS (**default**)
4. Study Creation
 - └ **System** **validates** configuration
 - └ **System** creates study record
 - └ Study appears in studies **list** with "**Ongoing**" status
 - └ **System** calculates estimated **end** date
5. Study Monitoring
 - └ Study progress tracked automatically
 - └ Hourly **data** quality monitoring
 - └ Real-time status updates

3. Review Holter Quality Data Workflow

Trigger: Clinician needs to assess study data quality

Flow Steps:

1. Access Study Progress
 - └─ Navigate **to** Studies section
 - └─ Locate target study in **list**
 - └─ Click "**Study Progress**" link
2. Study Progress Modal
 - └─ View overall study information
 - └─ Navigate **to** specific date using calendar
 - └─ Select date **to** view hourly breakdown
 - └─ Review hourly **data** quality
3. Quality Assessment
 - └─ Review "**Usable ECG / Total ECG**" metrics
 - └─ Interpret color-coded quality indicators
 - └─ Gray: 0% usable **data**
 - └─ Light Green: 1-24% usable **data**
 - └─ Green: 25-50% usable **data**
 - └─ Dark Green: >50% usable **data**
 - └─ Identify problematic time periods
 - └─ **Note** overall study quality trends
4. Quality-Based Actions
 - └─ **If** quality is poor:
 - └─ **Contact** patient **for** device adjustment
 - └─ Schedule device replacement
 - └─ Consider study extension
 - └─ **If** quality is **good**:
 - └─ **Continue** monitoring until study completion
5. Documentation
 - └─ Add clinical **notes** about quality assessment

4. Extend Study Workflow

Trigger: Study needs extension due to insufficient data or clinical requirements

Flow Steps:

1. Study Assessment
 - └─ Review current study progress
 - └─ Evaluate **data** quality
 - └─ Determine extension requirements
2. Study Modification
 - └─ Access study configuration (implied from UI)
 - └─ Modify **end** date
 - └─ Update study duration
 - └─ **Notify** patient of extension
3. **System** Updates
 - └─ Update estimated **end** date
 - └─ **Continue** **data** collection
 - └─ Reset quality monitoring **for** extended period

Note: Specific UI **for** study extension **not** visible in screenshots, but workflow inferred from study management patterns

5. Review ECG Strips per Symptom Workflow

Trigger: Patient reports symptoms that require ECG correlation

Flow Steps:

1. Access Symptom Logs
 - └─ Navigate to Symptom Log section
 - └─ Review list of reported symptoms
 - └─ Filter by status **if** needed
 - └─ Identify symptoms requiring review
2. Symptom Selection
 - └─ Click on specific symptom entry
 - └─ Note symptom details:
 - └─ Date **and** time
 - └─ Symptoms experienced
 - └─ Patient description
 - └─ Current status
 - └─ Click patient name to enter patient context
3. ECG Viewer Access
 - └─ Navigate to ECG VIEWER tab
 - └─ System automatically loads ECG data **for** symptom timepoint
 - └─ ECG viewer displays relevant time period
4. ECG Analysis
 - └─ Review ECG waveforms at symptom time
 - └─ Use filter controls to optimize display:
 - └─ Adjust frequency filtering
 - └─ Modify amplitude scaling
 - └─ Apply noise filtering
 - └─ Select appropriate channels
 - └─ Navigate through time periods:
 - └─ Use Previous/Next 4m controls
 - └─ Review expanded 10-second views
 - └─ Correlate with symptom timing
 - └─ Review contextual metrics:
 - └─ Blood pressure
 - └─ Activity level
 - └─ Other physiological data
 - └─ Temporal correlation
5. Clinical Annotation
 - └─ Click **"Add annotations" for** the symptom
 - └─ Complete observation checklist:
 - └─ Select relevant ECG findings
 - └─ Note arrhythmias **or** abnormalities
 - └─ Document technical observations
 - └─ Add clinical notes (1000 **char** limit)
 - └─ Provide interpretation (**for** reviewing physicians)
 - └─ Save annotations
6. Status Update
 - └─ Update symptom status:
 - └─ **"Pending review" → "MD to review"**
 - └─ **"MD to review" → "Reviewed by MD"**
 - └─ **"Reviewed by tech" (for** technical review)
 - └─ System tracks review workflow
7. Clinical Documentation
 - └─ Navigate to Notes section
 - └─ Create clinical note documenting findings
 - └─ Include symptom correlation
 - └─ Document clinical significance

6. Clinical Review and Diagnosis Workflow

Trigger: MD needs to review and provide clinical interpretation

Flow Steps:

1. Review Queue Access
 - |— Navigate **to** Symptom Log
 - |— Filter by **"MD to review"** status
 - |— Prioritize by clinical urgency
2. Symptom Review
 - |— Select symptom requiring MD review
 - |— Review technical annotations
 - |— Access ECG viewer **for** detailed analysis
 - |— Review patient **context** **and** history
3. Clinical Analysis
 - |— Correlate ECG findings with symptoms
 - |— Consider patient medical history
 - |— Evaluate clinical significance
 - |— Determine diagnosis **or** recommendations
4. Clinical Documentation
 - |— Update **interpretation** section in annotations
 - |— Provide clinical diagnosis
 - |— Document treatment recommendations
 - |— **Note** follow-up requirements
5. Status Completion
 - |— Update status **to** **"Reviewed by MD"**
 - |— **System** tracks completion
 - |— Triggers any required **noti**fications
6. Patient Communication
 - |— Generate patient report (implied)
 - |— Schedule follow-up **if** needed
 - |— Coordinate care team communication

7. Organization Management Workflow

Trigger: Need to manage team members and access

Flow Steps:

1. Access **Organization** Section
 - └─ Navigate **to** **Organization** in sidebar
2. Member Management
 - └─ Review current team members
 - └─ Check member roles **and** status
 - └─ Identify access needs
3. Add **New** Members
 - └─ Click **"Add Member"** button
 - └─ Enter member email address
 - └─ Assign appropriate role:
 - └─ Desk Admin
 - └─ MD
 - └─ Cardiac Tech
 - └─ **Organization** Admin
 - └─ Send invitation
 - └─ Track invitation status
4. Member Administration
 - └─ Monitor invitation acceptance
 - └─ Update member roles as needed
 - └─ Manage access permissions
 - └─ *Remove members when necessary*

Workflow Integration Points

Cross-Section Data Flow

Patient Creation → Study Setup → Data Collection → Symptom Logging → Clinical Review → Documentation

Each step creates data dependencies:

- ```

└─ Patient Profile enables Study Creation
└─ Active Study enables ECG Data Collection
└─ ECG Data enables Symptom Correlation
└─ Symptom Reports trigger Clinical Review
└─ Clinical Review generates Documentation

```

## Status Progression Workflows

### Symptom Status Progression:

Patient Reports → "Pending review" → "MD to review" → "Reviewed by MD"  
↘ "Reviewed by tech" ↗

### Study Status Progression:

```
Setup → "Ongoing" → "Completed" → "Closed by Staff"
 ↘ "Cancelled" (if needed)
```

### Patient Status Progression:

Registration → "Pending" → "Verified" → "Active" → "Archived"

## Quality Assurance Workflows

Data Quality Monitoring:  
Continuous ECG Collection → Hourly Quality Assessment → Quality Alerts → Intervention Actions

Clinical Quality Assurance:  
Technical Review → MD Review → Documentation → Quality Metrics → Process Improvement

## Error Handling and Edge Cases

### Common Error Scenarios

1. **Incomplete Patient Information:** Form validation prevents submission
2. **Study Setup Conflicts:** Date/time validation prevents scheduling conflicts
3. **Poor ECG Quality:** Quality monitoring triggers intervention workflows
4. **Missing Symptom Correlation:** System prompts for ECG review
5. **Incomplete Clinical Review:** Status tracking ensures completion

### Recovery Workflows

1. **Data Loss Recovery:** System maintains audit trails for data reconstruction
2. **Device Malfunction:** Study extension and device replacement workflows
3. **Clinical Escalation:** Urgent findings trigger immediate notification workflows
4. **System Downtime:** Offline data collection with sync capabilities

## Data Model and Backend Structure

---

### Database Schema Inference

Based on the UI elements and workflows observed, the following database structure can be inferred:

#### 1. Users Table

**Purpose:** Manages clinician accounts and authentication

```

CREATE TABLE users (
 id BIGINT PRIMARY KEY AUTO_INCREMENT,
 email VARCHAR(255) UNIQUE NOT NULL,
 password_hash VARCHAR(255) NOT NULL,
 first_name VARCHAR(100) NOT NULL,
 last_name VARCHAR(100) NOT NULL,
 role ENUM('desk_admin', 'md', 'cardiac_tech', 'organization_admin') NOT
 NULL,
 organization_id BIGINT NOT NULL,
 is_active BOOLEAN DEFAULT TRUE,
 mfa_enabled BOOLEAN DEFAULT FALSE,
 mfa_secret VARCHAR(255),
 created_at TIMESTAMP DEFAULT CURRENT_TIMESTAMP,
 updated_at TIMESTAMP DEFAULT CURRENT_TIMESTAMP ON UPDATE CURRENT_TIMESTAMP,
 last_login_at TIMESTAMP,

 FOREIGN KEY (organization_id) REFERENCES organizations(id)
);

```

**Key Fields Observed:** - Email: Used for login (simon.yang.ch@gmail.com) - Name: Displayed in header (Simon Yang) - Role: Determines access permissions (MD, Desk Admin, Cardiac Tech) - Organization: Links to clinic (Staging Clinic, ID: 373535) - MFA: Multi-factor authentication toggle visible

## 2. Organizations Table

**Purpose:** Manages clinic/organization information

```

CREATE TABLE organizations (
 id BIGINT PRIMARY KEY AUTO_INCREMENT,
 name VARCHAR(255) NOT NULL,
 organization_code VARCHAR(50) UNIQUE NOT NULL,
 address TEXT,
 phone VARCHAR(50),
 email VARCHAR(255),
 country VARCHAR(100),
 timezone VARCHAR(100),
 created_at TIMESTAMP DEFAULT CURRENT_TIMESTAMP,
 updated_at TIMESTAMP DEFAULT CURRENT_TIMESTAMP ON UPDATE CURRENT_TIMESTAMP
);

```

**Key Fields Observed:** - Name: "Staging Clinic" - Organization Code: "373535" - Multiple organizations supported (member invitations across organizations)

## 3. Patients Table

**Purpose:** Stores patient demographic and medical information

```

CREATE TABLE patients (
 id BIGINT PRIMARY KEY AUTO_INCREMENT,
 scla_account_id VARCHAR(255) UNIQUE,
 organization_id BIGINT NOT NULL,
 first_name VARCHAR(100) NOT NULL,
 last_name VARCHAR(100) NOT NULL,
 date_of_birth DATE NOT NULL,
 sex ENUM('male', 'female', 'unspecified'),
 email VARCHAR(255) NOT NULL,
 phone VARCHAR(50) NOT NULL,
 country VARCHAR(100) NOT NULL,

 -- Insurance Information
 health_insurance_policy VARCHAR(255) NOT NULL,
 insurance_number VARCHAR(100) NOT NULL,
 insurance_version_code VARCHAR(50),
 insurance_expiry_date DATE,

 -- Medical Information
 medical_record_number VARCHAR(100),
 assigned_physician_id BIGINT,

 -- Address Information
 street_address VARCHAR(255),
 city VARCHAR(100),
 province_state VARCHAR(100),
 postal_code VARCHAR(20),

 -- System Fields
 status ENUM('pending', 'verified', 'active', 'archived') DEFAULT 'pending',
 created_at TIMESTAMP DEFAULT CURRENT_TIMESTAMP,
 updated_at TIMESTAMP DEFAULT CURRENT_TIMESTAMP ON UPDATE CURRENT_TIMESTAMP,
 created_by BIGINT NOT NULL,

 FOREIGN KEY (organization_id) REFERENCES organizations(id),
 FOREIGN KEY (assigned_physician_id) REFERENCES users(id),
 FOREIGN KEY (created_by) REFERENCES users(id)
);

```

**Key Fields Observed:** - Personal: Yang, Simon, 1989-07-10, Male - Contact: simon.yang.ch@gmail.com, +41 78 795 00 09 - Insurance: Avenir (Groupe Mutuel), 7561234567897 - Address: Ackersteinstrasse 76, Zurich, 8049 - Status: "Verified" badge visible - Physician: Peter Wood (assigned)

#### 4. Studies Table

**Purpose:** Manages Holter monitoring studies

```

CREATE TABLE studies (
 id BIGINT PRIMARY KEY AUTO_INCREMENT,
 patient_id BIGINT NOT NULL,
 organization_id BIGINT NOT NULL,
 study_type ENUM('holter') DEFAULT 'holter',
 duration_type ENUM('24_hour', '48_hour', '72_hour', '7_day', '14_day') NOT
NULL,

 -- Study Timeline
 start_date DATE NOT NULL,
 start_time TIME NOT NULL,
 estimated_end_date DATE NOT NULL,
 actual_end_date DATE,

 -- Study Configuration
 setup_type VARCHAR(100) DEFAULT 'in_clinic_setup',
 garment_type VARCHAR(100) DEFAULT 'skipp_band',
 garment_size VARCHAR(20) DEFAULT 'xs',

 -- Study Status
 status ENUM('ongoing', 'completed', 'cancelled', 'closed_by_staff') DEFAULT
'ongoing',

 -- Quality Metrics
 total_ecg_minutes INT DEFAULT 0,
 usable_ecg_minutes INT DEFAULT 0,
 average_daily_usable_hours DECIMAL(4,2) DEFAULT 0.00,

 -- Associated Data
 pod_id VARCHAR(255),
 edf_file_path VARCHAR(500),

 -- System Fields
 created_at TIMESTAMP DEFAULT CURRENT_TIMESTAMP,
 updated_at TIMESTAMP DEFAULT CURRENT_TIMESTAMP ON UPDATE CURRENT_TIMESTAMP,
 created_by BIGINT NOT NULL,

 FOREIGN KEY (patient_id) REFERENCES patients(id),
 FOREIGN KEY (organization_id) REFERENCES organizations(id),
 FOREIGN KEY (created_by) REFERENCES users(id)
);

```

**Key Fields Observed:** - Duration: 14-days, 48-hours, 24-hours options - Dates: 2025-06-24 10:21 start, 2025-07-08 estimated end - Status: "Ongoing", "Closed by Staff", "Cancelled" - Quality: "Avg. 0 hrs of usable ECG/day" - Pod Association: "31067601890"

## 5. Study\_Quality Table

**Purpose:** Stores hourly ECG data quality metrics

```

CREATE TABLE study_quality (
 id BIGINT PRIMARY KEY AUTO_INCREMENT,
 study_id BIGINT NOT NULL,
 date DATE NOT NULL,
 hour TINYINT NOT NULL, -- 0-23

 -- Quality Metrics
 total_minutes TINYINT DEFAULT 0, -- 0-60
 usable_minutes TINYINT DEFAULT 0, -- 0-60
 quality_percentage DECIMAL(5,2) DEFAULT 0.00,
 quality_category ENUM('none', 'low', 'medium', 'high') DEFAULT 'none',

 -- Technical Metrics
 signal_quality_score DECIMAL(5,2),
 noise_level DECIMAL(5,2),
 artifact_count INT DEFAULT 0,

 -- System Fields
 created_at TIMESTAMP DEFAULT CURRENT_TIMESTAMP,
 updated_at TIMESTAMP DEFAULT CURRENT_TIMESTAMP ON UPDATE CURRENT_TIMESTAMP,

 FOREIGN KEY (study_id) REFERENCES studies(id),
 UNIQUE KEY unique_study_hour (study_id, date, hour)
);

```

**Key Fields Observed:** - Hourly breakdown: 00:00-23:00 time slots - Quality metrics: "0 mins / 0 mins" (usable/total) - Quality categories: Color-coded (0%, 1-24%, 25-50%, >50%) - Date-specific tracking: Calendar navigation for specific dates

## 6. Symptoms Table

**Purpose:** Stores patient-reported symptoms

```

CREATE TABLE symptoms (
 id BIGINT PRIMARY KEY AUTO_INCREMENT,
 patient_id BIGINT NOT NULL,
 study_id BIGINT,

 -- Symptom Information
 symptom_date TIMESTAMP NOT NULL,
 reported_date TIMESTAMP NOT NULL,
 symptoms_experienced TEXT NOT NULL,
 severity TINYINT, -- 1-10 scale
 duration VARCHAR(100),
 triggers TEXT,
 patient_notes TEXT,

 -- Clinical Review
 status ENUM('pending_review', 'md_to_review', 'reviewed_by_md',
'reviewed_by_tech') DEFAULT 'pending_review',
 assigned_physician_id BIGINT,

 -- ECG Correlation
 ecg_start_time TIMESTAMP,
 ecg_end_time TIMESTAMP,
 has_routine_ecg BOOLEAN DEFAULT FALSE,

 -- System Fields
 created_at TIMESTAMP DEFAULT CURRENT_TIMESTAMP,
 updated_at TIMESTAMP DEFAULT CURRENT_TIMESTAMP ON UPDATE CURRENT_TIMESTAMP,

 FOREIGN KEY (patient_id) REFERENCES patients(id),
 FOREIGN KEY (study_id) REFERENCES studies(id),
 FOREIGN KEY (assigned_physician_id) REFERENCES users(id)
);

```

**Key Fields Observed:** - Timing: "2025-06-24 10:38" symptom and reported dates - Symptoms: "Sweating", "testing", "study started", "Routine ECG" - Status: "Pending review" with dropdown options - Severity: "5/10" scale - Duration: "On Going" - Triggers: "Caffeine"

## 7. ECG\_Data Table

**Purpose:** Links to ECG binary data stored in S3



```

CREATE TABLE ecg_data (
 id BIGINT PRIMARY KEY AUTO_INCREMENT,
 study_id BIGINT NOT NULL,
 patient_id BIGINT NOT NULL,

 -- Temporal Information
 start_timestamp TIMESTAMP NOT NULL,
 end_timestamp TIMESTAMP NOT NULL,
 duration_seconds INT NOT NULL,

 -- ECG Channels
 channel_1_s3_path VARCHAR(500),
 channel_2_s3_path VARCHAR(500),
 channel_3_s3_path VARCHAR(500),

 -- Data Quality
 signal_quality ENUM('poor', 'fair', 'good', 'excellent'),
 noise_level DECIMAL(5,2),
 artifact_flags JSON,

 -- Processing Information
 sampling_rate INT DEFAULT 250, -- Hz
 amplitude_scale DECIMAL(8,4) DEFAULT 10.0, -- mm/mV
 frequency_filter_low DECIMAL(6,3) DEFAULT 0.05, -- Hz
 frequency_filter_high DECIMAL(6,3) DEFAULT 45.0, -- Hz

 -- System Fields
 created_at TIMESTAMP DEFAULT CURRENT_TIMESTAMP,
 processed_at TIMESTAMP,

 FOREIGN KEY (study_id) REFERENCES studies(id),
 FOREIGN KEY (patient_id) REFERENCES patients(id),

 INDEX idx_study_timestamp (study_id, start_timestamp),
 INDEX idx_patient_timestamp (patient_id, start_timestamp)
);

```

**Key Fields Observed:** - Time ranges: "10:34:30 - 10:38:30" (4-minute segments) - Channels: 1, 2, 3 with radio button selection - Filters: "0.05-45 Hz (Standard)", "Off" noise filter, "10mm/mV" amplitude - S3 storage: Binary ECG data stored externally

## 8. Symptom\_Annotations Table

**Purpose:** Stores clinical annotations for symptoms

```

CREATE TABLE symptom_annotations (
 id BIGINT PRIMARY KEY AUTO_INCREMENT,
 symptom_id BIGINT NOT NULL,
 annotated_by BIGINT NOT NULL,

 -- Observation Checkboxes (stored as JSON for flexibility)
 observations JSON, -- Array of selected observations

 -- Text Fields
 clinical_notes TEXT, -- 1000 character limit
 interpretation TEXT, -- 1000 character limit (for reviewing physicians only)

 -- System Fields
 created_at TIMESTAMP DEFAULT CURRENT_TIMESTAMP,
 updated_at TIMESTAMP DEFAULT CURRENT_TIMESTAMP ON UPDATE CURRENT_TIMESTAMP,

 FOREIGN KEY (symptom_id) REFERENCES symptoms(id),
 FOREIGN KEY (annotated_by) REFERENCES users(id)
);

```

**Key Fields Observed:** - Observations: Checkboxes for 1AVB, 2AVB1, 2AVB2, 3AVB, AF, Artifact, AFL, BBB, Diagnosis, EAT, IRD, IVR, JR, MAT, No ECG, NNS, NNV, NS, NV, Pause, SA, SR, SVT, VFib, VT, WAP, Other - Notes: 1000 character limit - Interpretation: 1000 character limit (physicians only)

## 9. Clinical\_Notes Table

**Purpose:** Stores clinical documentation

```

CREATE TABLE clinical_notes (
 id BIGINT PRIMARY KEY AUTO_INCREMENT,
 patient_id BIGINT NOT NULL,
 created_by BIGINT NOT NULL,

 -- Note Content
 title VARCHAR(255) NOT NULL,
 content TEXT NOT NULL,
 note_type ENUM('clinical', 'administrative', 'technical') DEFAULT 'clinical',

 -- Rich Text Formatting (stored as HTML or markdown)
 formatted_content TEXT,

 -- System Fields
 created_at TIMESTAMP DEFAULT CURRENT_TIMESTAMP,
 -- Notes cannot be edited or deleted once created (per UI notice)

 FOREIGN KEY (patient_id) REFERENCES patients(id),
 FOREIGN KEY (created_by) REFERENCES users(id)
);

```

**Key Fields Observed:** - Rich text editor with formatting options - Title field - Immutable once created (per UI warning) - Date/time auto-populated

## 10. Contextual\_Metrics Table

**Purpose:** Stores physiological data for ECG correlation

```
CREATE TABLE contextual_metrics (
 id BIGINT PRIMARY KEY AUTO_INCREMENT,
 patient_id BIGINT NOT NULL,
 study_id BIGINT,

 -- Temporal Information
 measurement_timestamp TIMESTAMP NOT NULL,

 -- Physiological Metrics
 blood_pressure_systolic INT,
 blood_pressure_diastolic INT,
 oxygen_saturation DECIMAL(5,2),
 body_weight DECIMAL(6,2),
 glucose_level DECIMAL(6,2),
 body_temperature DECIMAL(4,2),
 activity_steps INT,
 last_posture VARCHAR(50),

 -- Data Source
 source_device VARCHAR(100),
 measurement_type ENUM('manual', 'automatic', 'imported'),

 -- System Fields
 created_at TIMESTAMP DEFAULT CURRENT_TIMESTAMP,

 FOREIGN KEY (patient_id) REFERENCES patients(id),
 FOREIGN KEY (study_id) REFERENCES studies(id),

 INDEX idx_patient_timestamp (patient_id, measurement_timestamp)
);
```

**Key Fields Observed:** - Blood pressure: "mmHg (closest before/after)" - Oxygen saturation: "mmHg (closest before/after)" - Body weight: "kg (closest before/after)" - Glucose: "mmol/L (closest before/after)" - Temperature: "°C" - Activity: "steps", "last posture"

## Backend Architecture Inference

### API Structure

Based on the UI interactions, the following API endpoints can be inferred:

#### Authentication & Users:

POST /api/auth/login  
POST /api/auth/logout  
GET /api/auth/me  
PUT /api/auth/mfa/enable  
PUT /api/auth/mfa/disable

#### Organizations:

GET /api/organizations/{id}/members  
POST /api/organizations/{id}/members/invite  
DELETE /api/organizations/{id}/members/{userId}

#### Patients:

GET /api/patients  
POST /api/patients  
GET /api/patients/{id}  
PUT /api/patients/{id}  
DELETE /api/patients/{id}/archive  
POST /api/patients/{id}/scla-login-code

#### Studies:

GET /api/studies  
POST /api/studies  
GET /api/studies/{id}  
PUT /api/studies/{id}  
GET /api/studies/{id}/progress  
GET /api/studies/{id}/quality/{date}  
GET /api/studies/{id}/edf-download

#### Symptoms:

GET /api/symptoms  
GET /api/symptoms/{id}  
PUT /api/symptoms/{id}/status  
GET /api/symptoms/{id}/history  
POST /api/symptoms/{id}/annotations

#### ECG Data:

GET /api/ecg/{studyId}/data  
GET /api/ecg/{studyId}/channels/{channel}  
GET /api/ecg/{studyId}/timerange/{start}/{end}  
POST /api/ecg/filters/apply

#### Notes:

GET /api/notes  
POST /api/notes  
GET /api/notes/{id}

#### Contextual Metrics:

GET /api/metrics/{patientId}/timerange/{start}/{end}  
POST /api/metrics/{patientId}

## Data Storage Architecture

Primary Database (MySQL/PostgreSQL):

- └─ User management **and** authentication
- └─ Patient demographics **and** medical records
- └─ Study configuration **and** metadata
- └─ Symptom logs **and** annotations
- └─ Clinical notes **and** documentation
- └─ Quality metrics **and** aggregations

S3 **Object** Storage:

- └─ Raw ECG binary data files
- └─ Processed ECG segments
- └─ EDF **export** files
- └─ System backups

Redis Cache:

- └─ Session management
- └─ Real-time ECG data cache
- └─ Quality metric aggregations
- └─ API response caching

Time-Series Database (InfluxDB/TimescaleDB):

- └─ High-frequency ECG data points
- └─ Contextual physiological metrics
- └─ Device telemetry data
- └─ Quality monitoring metrics

## Integration Points

External Systems:

- └─ SCLA Mobile App (patient-facing)
- └─ ECG Device APIs (data collection)
- └─ Healthcare Information Systems (HL7/FHIR)
- └─ Insurance Verification Services
- └─ Notification Services (email/SMS)

Internal Services:

- └─ ECG Processing Engine
- └─ Quality Assessment Service
- └─ Clinical Decision Support
- └─ Report Generation Service
- └─ Audit Logging Service

# Data Relationships and Constraints

## Primary Relationships

```
Organizations (1) → (N) Users
Organizations (1) → (N) Patients
Patients (1) → (N) Studies
Studies (1) → (N) ECG_Data
Studies (1) → (N) Study_Quality
Patients (1) → (N) Symptoms
Symptoms (1) → (N) Symptom_Annotations
Patients (1) → (N) Clinical_Notes
Patients (1) → (N) Contextual_Metrics
```

## Business Rules and Constraints

1. **Patient-Study Relationship:** One patient can have multiple studies, but each study belongs to one patient
2. **Study Quality:** Quality data is generated hourly for active studies
3. **Symptom Correlation:** Symptoms must be correlated with ECG data within study timeframe
4. **Clinical Review:** Only MDs can provide clinical interpretations
5. **Data Immutability:** Clinical notes cannot be modified once created
6. **Organization Isolation:** Users can only access data within their organization
7. **Study Lifecycle:** Studies progress through defined status states
8. **Quality Thresholds:** Minimum quality requirements for clinical validity

# Implementation Guidance

---

## Technology Stack Recommendations

### Frontend Technologies

```
Framework: React.js with TypeScript
├─ UI Library: Material-UI or Ant Design
├─ State Management: Redux Toolkit or Zustand
├─ Routing: React Router v6
├─ Forms: React Hook Form with Yup validation
├─ Charts: Chart.js or D3.js for ECG visualization
├─ Rich Text: Draft.js or TinyMCE for notes editor
└─ Date/Time: date-fns or moment.js
```

#### ECG Visualization:

```
├─ Canvas-based rendering for performance
├─ WebGL for real-time data processing
├─ Custom waveform components
└─ Multi-channel synchronization
```

### Backend Technologies

```
Framework: Node.js with Express.js or Python with FastAPI
├─ Database: PostgreSQL with TimescaleDB extension
├─ ORM: Prisma (Node.js) or SQLAlchemy (Python)
├─ Authentication: JWT with refresh tokens
├─ File Storage: AWS S3 or compatible object storage
├─ Caching: Redis for session and data caching
├─ Queue: Bull (Node.js) or Celery (Python) for background jobs
└─ Real-time: Socket.io or WebSockets for live updates
```

#### ECG Processing:

```
├─ Signal processing libraries (scipy, numpy)
├─ Real-time data streaming
├─ Quality assessment algorithms
└─ Artifact detection and filtering
```

## Infrastructure

### Deployment:

- Containerization: Docker with docker-compose
- Orchestration: Kubernetes or Docker Swarm
- Load Balancing: NGINX or AWS ALB
- Monitoring: Prometheus + Grafana
- Logging: ELK Stack (Elasticsearch, Logstash, Kibana)
- CI/CD: GitHub Actions or GitLab CI

### Security:

- HTTPS/TLS encryption
- HIPAA compliance measures
- Data encryption at rest and in transit
- Regular security audits
- Access logging and monitoring

## Component Implementation Details

### 1. ECG Viewer Component

#### Core Requirements:

```
interface ECGViewerProps {
 studyId: string;
 patientId: string;
 timeRange: {
 start: Date;
 end: Date;
 };
 channels: number[];
 filters: ECGFilters;
 onAnnotationAdd: (annotation: Annotation) => void;
}

interface ECGFilters {
 frequencyLow: number; // 0.05 Hz default
 frequencyHigh: number; // 45 Hz default
 mainsNoiseFilter: boolean;
 amplitude: number; // 10 mm/mV default
}
```

**Implementation Considerations:** - Canvas-based rendering for smooth scrolling and zooming - Efficient data loading with pagination for large datasets - Real-time filter application without server round-trips - Multi-channel synchronization with precise time alignment - Annotation overlay system with click-to-annotate functionality

### 2. Data Table Component

#### Reusable Table Structure:



```

interface DataTableProps<T> {
 data: T[];
 columns: ColumnDefinition<T>[];
 pagination: PaginationConfig;
 sorting: SortConfig;
 filtering: FilterConfig;
 actions?: ActionConfig<T>[];
}

interface ColumnDefinition<T> {
 key: keyof T;
 title: string;
 sortable: boolean;
 filterable: boolean;
 render?: (value: any, record: T) => React.ReactNode;
}

```

**Features to Implement:** - Server-side pagination for large datasets - Multi-column sorting with visual indicators - Advanced filtering with multiple criteria - Inline editing for status changes - Bulk actions for multiple selections

### 3. Modal System

#### Modal Management:

```

interface ModalConfig {
 id: string;
 component: React.ComponentType<any>;
 props: any;
 size: 'small' | 'medium' | 'large' | 'fullscreen';
 closable: boolean;
 maskClosable: boolean;
}

// Global modal state management
const useModalStore = () => {
 const [modals, setModals] = useState<ModalConfig[]>([]);

 const openModal = (config: ModalConfig) => { /* ... */ };
 const closeModal = (id: string) => { /* ... */ };
 const closeAllModals = () => { /* ... */ };

 return { modals, openModal, closeModal, closeAllModals };
};

```

### 4. Form Validation System

#### Validation Schema Examples:

```

// Patient Creation Form
const patientValidationSchema = yup.object({
 firstName: yup.string().required('First name is required'),
 lastName: yup.string().required('Last name is required'),
 dateOfBirth: yup.date().required('Date of birth is required'),
 email: yup.string().email('Invalid email').required('Email is required'),
 phone: yup.string().required('Phone number is required'),
 healthInsurancePolicy: yup.string().required('Insurance policy is required'),
 insuranceNumber: yup.string().required('Insurance number is required'),
});

// Study Creation Form
const studyValidationSchema = yup.object({
 duration: yup.string().oneOf(['24_hour', '48_hour', '72_hour', '7_day', '14_day']).required(),
 startDate: yup.date().min(new Date(), 'Start date cannot be in the past').required(),
 startTime: yup.string().matches(/^([0-1]?[0-9]|2[0-3]):[0-5][0-9]$/, 'Invalid time format').required(),
});

```

## API Implementation Specifications

### Authentication & Authorization

#### JWT Token Structure:

```

interface JWTPayload {
 userId: string;
 organizationId: string;
 role: UserRole;
 permissions: string[];
 iat: number;
 exp: number;
}

// Role-based access control
const permissions = {
 desk_admin: ['patients:read', 'patients:write', 'studies:read'],
 md: ['patients:read', 'studies:read', 'symptoms:review', 'notes:write'],
 cardiac_tech: ['studies:read', 'ecg:analyze', 'symptoms:review'],
 organization_admin: ['users:manage', 'organization:manage']
};

```

### Real-time Data Handling

#### WebSocket Events:

```

// Client-side event handling
interface ECGDataEvent {
 studyId: string;
 timestamp: Date;
 channels: {
 channel1: number[];
 channel2: number[];
 channel3: number[];
 };
 quality: QualityMetrics;
}

// Server-side event emission
socket.emit('ecg:data', {
 studyId: study.id,
 timestamp: new Date(),
 channels: processedECGData,
 quality: qualityAssessment
});

```

## Data Processing Pipeline

### ECG Quality Assessment:

```

def assess_ecg_quality(ecg_data, sampling_rate=250):
 """
 Assess ECG signal quality based on multiple criteria
 """
 quality_score = 0

 # Signal-to-noise ratio
 snr = calculate_snr(ecg_data)
 quality_score += min(snr / 20, 0.3) # Max 30% from SNR

 # Artifact detection
 artifact_percentage = detect_artifacts(ecg_data)
 quality_score += max(0, 0.3 - artifact_percentage) # Max 30% from artifacts

 # Heart rate variability
 hrv_score = assess_hrv_quality(ecg_data, sampling_rate)
 quality_score += hrv_score * 0.2 # Max 20% from HRV

 # Baseline stability
 baseline_score = assess_baseline_stability(ecg_data)
 quality_score += baseline_score * 0.2 # Max 20% from baseline

 return min(quality_score, 1.0)

```

# Security Implementation

## Data Encryption

### Encryption Strategy:

```
// Patient data encryption
const encryptPatientData = (data: PatientData): EncryptedPatientData => {
 const sensitiveFields = ['ssn', 'medicalRecordNumber', 'insuranceNumber'];
 const encrypted = { ...data };

 sensitiveFields.forEach(field => {
 if (encrypted[field]) {
 encrypted[field] = encrypt(encrypted[field], process.env.ENCRIPTION_KEY);
 }
 });

 return encrypted;
};

// ECG data integrity
const generateECGChecksum = (ecgData: ECGData): string => {
 return crypto
 .createHash('sha256')
 .update(JSON.stringify(ecgData))
 .digest('hex');
};
```

## Audit Logging

### Audit Trail Implementation:

```
interface AuditLog {
 id: string;
 userId: string;
 action: string;
 resourceType: string;
 resourceId: string;
 changes: any;
 ipAddress: string;
 userAgent: string;
 timestamp: Date;
}

const logAuditEvent = async (event: Omit<AuditLog, 'id' | 'timestamp'>) => {
 await auditLogRepository.create({
 ...event,
 id: generateUUID(),
 timestamp: new Date()
 });
};
```

# Performance Optimization

## ECG Data Optimization

### Data Compression and Streaming:

```
// Client-side ECG data management
class ECGDataManager {
 private cache = new Map<string, ECGSegment>();
 private maxCacheSize = 100; // segments

 async loadECGSegment(studyId: string, startTime: Date, duration: number):
 Promise<ECGSegment> {
 const cacheKey = `${studyId}-${startTime.getTime()}-${duration}`;

 if (this.cache.has(cacheKey)) {
 return this.cache.get(cacheKey)!;
 }

 const segment = await this.fetchECGData(studyId, startTime, duration);
 this.addToCache(cacheKey, segment);

 return segment;
 }

 private addToCache(key: string, segment: ECGSegment) {
 if (this.cache.size >= this.maxCacheSize) {
 const firstKey = this.cache.keys().next().value;
 this.cache.delete(firstKey);
 }
 this.cache.set(key, segment);
 }
}
```

## Database Optimization

### Query Optimization Strategies:

```
-- Optimized symptom query with proper indexing
CREATE INDEX idx_symptoms_patient_date ON symptoms(patient_id, symptom_date
DESC);
CREATE INDEX idx_symptoms_status_date ON symptoms(status, symptom_date DESC);

-- Optimized ECG data query for time-range searches
CREATE INDEX idx_ecg_study_time ON ecg_data(study_id, start_timestamp,
end_timestamp);

-- Partitioned study_quality table for better performance
CREATE TABLE study_quality_2025 PARTITION OF study_quality
FOR VALUES FROM ('2025-01-01') TO ('2026-01-01');
```

# Testing Strategy

## Unit Testing

### Component Testing Examples:

```
// ECG Viewer component test
describe('ECGViewer', () => {
 it('should render ECG waveforms correctly', async () => {
 const mockData = generateMockECGData();
 render(<ECGViewer studyId="123" data={mockData} />);

 expect(screen.getByTestId('ecg-canvas')).toBeInTheDocument();
 expect(screen.getByText('Channel 1')).toBeInTheDocument();
 });

 it('should apply filters correctly', async () => {
 const mockData = generateMockECGData();
 const { rerender } = render(<ECGViewer studyId="123" data={mockData} />);

 const filters = { frequencyLow: 0.1, frequencyHigh: 40, amplitude: 15 };
 rerender(<ECGViewer studyId="123" data={mockData} filters={filters} />);

 // Verify filter application
 expect(mockFilterFunction).toHaveBeenCalledWith(mockData, filters);
 });
});
```

## Integration Testing

### API Integration Tests:

```
describe('Patient API', () => {
 it('should create patient with valid data', async () => {
 const patientData = {
 firstName: 'John',
 lastName: 'Doe',
 dateOfBirth: '1990-01-01',
 email: 'john.doe@example.com',
 // ... other required fields
 };

 const response = await request(app)
 .post('/api/patients')
 .set('Authorization', `Bearer ${validToken}`)
 .send(patientData)
 .expect(201);

 expect(response.body.id).toBeDefined();
 expect(response.body.status).toBe('pending');
 });
});
```

# Deployment Configuration

## Docker Configuration

### Dockerfile Example:

```
Frontend Dockerfile
FROM node:18-alpine AS builder
WORKDIR /app
COPY package*.json ./
RUN npm ci --only=production
COPY . .
RUN npm run build

FROM nginx:alpine
COPY --from=builder /app/dist /usr/share/nginx/html
COPY nginx.conf /etc/nginx/nginx.conf
EXPOSE 80
CMD ["nginx", "-g", "daemon off;"]
```

### Docker Compose Configuration:

```
version: '3.8'
services:
 frontend:
 build: ./frontend
 ports:
 - "3000:80"
 depends_on:
 - backend

 backend:
 build: ./backend
 ports:
 - "8000:8000"
 environment:
 - DATABASE_URL=postgresql://user:pass@db:5432/mvcp
 - REDIS_URL=redis://redis:6379
 depends_on:
 - db
 - redis

 db:
 image: timescale/timescaledb:latest-pg14
 environment:
 - POSTGRES_DB=mvcp
 - POSTGRES_USER=user
 - POSTGRES_PASSWORD=pass
 volumes:
 - postgres_data:/var/lib/postgresql/data

 redis:
 image: redis:alpine
 volumes:
 - redis_data:/data

volumes:
 postgres_data:
 redis_data:
```

## Monitoring and Maintenance

### Health Checks

#### Application Health Monitoring:



```
// Health check endpoint
app.get('/health', async (req, res) => {
 const health = {
 status: 'ok',
 timestamp: new Date().toISOString(),
 services: {
 database: await checkDatabaseHealth(),
 redis: await checkRedisHealth(),
 s3: await checkS3Health(),
 ecgProcessor: await checkECGProcessorHealth()
 }
 };

 const isHealthy = Object.values(health.services).every(service =>
 service.status === 'ok');
 res.status(isHealthy ? 200 : 503).json(health);
});
```

## Performance Monitoring

### Key Metrics to Track:

```
const metrics = {
 // Application Performance
 responseTime: 'Average API response time',
 throughput: 'Requests per second',
 errorRate: 'Percentage of failed requests',

 // ECG Processing
 ecgProcessingLatency: 'Time to process ECG segments',
 qualityAssessmentTime: 'Time to assess data quality',
 dataIngestionRate: 'ECG data points per second',

 // User Experience
 pageLoadTime: 'Frontend page load times',
 ecgRenderTime: 'Time to render ECG waveforms',
 modalOpenTime: 'Time to open complex modals',

 // Business Metrics
 activeStudies: 'Number of ongoing studies',
 dailySymptomReports: 'Symptoms reported per day',
 clinicalReviewTime: 'Average time for clinical review'
};
```

This comprehensive documentation provides all the necessary information for an AI system to reconstruct the MVCP application from scratch, including detailed technical specifications, implementation guidance, and operational considerations.